**Pertussis
DNA Amplification Assay**

DNA Amplification Assay for the Detection of *Bordetella pertussis* in nasopharyngeal swab samples

REF 280750

IVD In vitro diagnostic medical device

INTENDED USE

The *illumigene* Pertussis DNA Amplification Assay, performed on the *illumipro-10™*, is a qualitative in vitro diagnostic test for the direct detection of *Bordetella pertussis* in human nasopharyngeal swab samples taken from patients suspected of having respiratory tract infection attributable to *Bordetella pertussis*.

The *illumigene* Pertussis assay utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect *Bordetella pertussis* by targeting the IS481 insertional element of the *Bordetella pertussis* genome. The IS481 insertional element can also be found in *Bordetella holmesii* and some *Bordetella bronchiseptica* strains. While respiratory illness caused by *B. bronchiseptica* is infrequent in humans, mixed outbreaks involving *B. holmesii* and *B. pertussis* have been reported. Positive results reported for DNA amplification assays targeting IS481 may be consistent with the presence of either *B. holmesii* or *B. pertussis*. When clinical factors suggest that *B. pertussis* may not be the cause of respiratory infection, other clinically appropriate investigation(s) should be carried out in accordance with published guidelines.

Results from the *illumigene* Pertussis assay should be used in conjunction with information obtained during the patient's clinical evaluation as an aid in diagnosis of respiratory infection.

illumigene Pertussis is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The *illumigene* Pertussis DNA Amplification Assay is based on loop-mediated amplification (LAMP)^{1,2} technology. The assay targets a 198 base pair (bp) sequence of the *Bordetella pertussis* genome residing in a region of the IS481 insertional element sequence.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of this amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian *illumipro-10* Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The *illumigene* Pertussis kit includes *illumigene* Pertussis Assay Control/Negative Control Reagent, *illumigene* Pertussis Test Devices, *illumigene* Pertussis Sample Buffer and Mineral Oil. The *illumigene* Assay Control/Negative Control, used for specimen dilution and preparation, is a Tris-buffered solution containing formalin-treated *E. coli* harboring *Staphylococcus aureus* DNA. The *illumigene* Pertussis Test Device contains one lyophilized amplification bead in each of two chambers: a TEST chamber with IS481-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA in the Assay Control/Negative Control Reagent and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is added to the Assay Control/Negative Control Reagent and combined with the *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient sample allows for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The Control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The *illumipro-10* monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_f). The *illumipro-10* calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Bordetella pertussis is a human pathogen exclusively responsible for the endemic respiratory disease whooping cough (pertussis); a publically reportable disease that impacted over 200,000 people worldwide in 2012. Whooping cough is currently the only disease in the US that displays an increasing trend of reportable cases and substantial morbidity despite the availability of a vaccine.³ The US Centers for Disease Control (CDC) reported approximately 48,300 cases of pertussis in 2012, an average of 15.4 cases per 100,000 people in the United States.⁴ Recent pertussis outbreaks have indicated waning immunity in children after receiving the final booster acellular vaccine dose, allowing for susceptibility of pertussis infection in both vaccinated and unvaccinated populations of children ages 7-11.^{5,6}

Pertussis outbreaks are particularly challenging to control as the early stages of disease resemble other respiratory infections and patients can remain highly infectious for up to 5 weeks after onset of symptoms.³ In addition, pertussis is not officially diagnosed until the presence of the characteristic "whooping" cough that is displayed almost 2 weeks post-symptom onset.³ Clinical diagnosis of pertussis can be performed using culture, serology, or nucleic acid amplification tests (NAAT). While culture is highly specific, sensitivity is low, results can take up to 7 days, and the viable cells required for culture decrease with disease progression (after 2 weeks of symptom onset).⁷ Unlike culture, NAAT assays can provide rapid test results and do not require viable bacteria. However, bacterial DNA must still be present in the nasopharynx (up to 4 weeks after cough onset) in order to prevent falsely-negative results.⁸ In addition, only symptomatic patients with cough should be tested by NAAT in order to prevent false positive results by asymptomatic close contacts.⁹ Serology can only be used for diagnosis in late stages of disease, approximately 2-8 weeks after cough onset.

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. ***illumigene* Pertussis Assay Control/Negative Control Reagent:** Tris-buffered solution containing formalin-treated *E. coli* harboring plasmid containing a segment of the *S. aureus* genome and sodium azide (0.09%) as a preservative.
2. ***illumigene* Pertussis Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates), and either IS481-specific primers (TEST Chamber) or control primers (CONTROL Chamber).
3. ***illumigene* Pertussis Sample Buffer:** Tris-EDTA solution containing sodium azide (0.09%) as a preservative.
4. **Mineral Oil**

MATERIALS PROVIDED SEPARATELY

illumigene Pertussis External Control Kit, Catalog Number: 279930

MATERIALS NOT PROVIDED

1. Disposable latex gloves, powder free
 2. DNase/RNase-free, aerosol resistant pipette tips
 3. Specimen Collection and transport system
- Nasopharyngeal Swabs:** Polyester (Minimum Capacity 18µL, e.g. Puritan Medical Products Catalog 25-801-D50); Flocked nylon (Minimum Capacity 31 µL, e.g. Puritan Medical Products Catalog 25-201-R50); or Rayon (Minimum Capacity 69 µL, e.g. Copan 503CS01).
- Non-nutritive Transport Medium:** Liquid Amies, without charcoal; or Liquid Stuart (Maximum Volume: 1.2 mL)

EQUIPMENT NOT PROVIDED

1. Dry-bath with 12 mm heat block capable of 95 C
2. Digital thermometer with Max/Min Temperature Memory (eg, Traceable® Lollipop™ Waterproof/Shockproof Thermometer)
3. Vortex mixer
4. Interval timer
5. Micropipette capable of dispensing 50 µL
6. *illumipro-10*, Meridian Bioscience, Inc. Catalog Number: 610172

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange Assay Control/Negative Control Reagent or Test Devices between lots. Sample Buffer and Mineral Oil are interchangeable provided they are within assigned expiration dates when used.
3. Follow Biosafety Level 2 and Good Laboratory practices during testing.⁹ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
4. Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
5. Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.¹⁰
6. The *illumigene* Pertussis Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
7. The *illumigene* Pertussis Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
8. Dispose of used *illumigene* Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

RISK AND SAFETY PHRASES

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-30 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-30 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION AND PREPARATION

Sample type: Nasopharyngeal swabs

Sample Collection: Nasopharyngeal swab specimen collection should be performed in accordance with institutional procedures for collection of clinical specimens for *Bordetella pertussis* infection. Nasopharyngeal swab samples should be collected with suitable swab types (e.g. Polyester, Flocked Nylon or Rayon).

Place swab(s) in non-nutritive transport medium (eg, Liquid Amies, without charcoal or Liquid Stuart) or store unpreserved in a sterile tube without medium.

Unpreserved samples or samples stored in transport media should be tested as soon as possible, but may be held at room temperature (21–30 C) for up to 5 days or refrigerated (2–8 C) for up to 7 days prior to testing. Do not freeze samples.

SPECIMEN PREPARATION:

NOTE: Ensure that the *illumipro-10* is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

1. Place the swab specimen in a labeled Sample Buffer tube. Cut the swab shaft to ensure sample fits in the tube. Elute the sample by vortexing for 45-60 seconds. Remove and discard the swab. **Eluted samples may be held at room temperature (21-30 C) for up to 48 hours or refrigerated (2-8 C) for up to 7 days prior to testing.**
2. Add 50 µL of eluted sample to a labeled Assay Control/Negative Control Tube and cap.
3. Repeat Specimen Preparation Steps for all samples to be processed.
4. Vortex each Assay Control/Negative Control Tube containing eluted sample for approximately 10 seconds.
5. Heat each Sample/Control mixture in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
6. Remove each Sample/Control tube from the dry-bath/heat block. Heat treated samples may be held at room temperature (21-30 C) for up to 15 minutes prior to testing.
7. Vortex for approximately 10 seconds.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single *illumipro-10* run.

1. Remove 1 *illumigene* Pertussis Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagent will not fall out upon opening.
2. Transfer 50 µL of the heat-treated sample to the TEST chamber (White Bead) of the *illumigene* Test Device. Take care not to introduce air to the reaction mixture. Using a new pipette tip, transfer 50 µL of the heat-treated sample to the CONTROL chamber (Yellow Bead) of the *illumigene* Test Device. Take care not to introduce air to the reaction mixture.
3. Add 1 drop of Mineral Oil to both the TEST chamber and CONTROL chamber. Close the *illumigene* Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the Test Device for rehydration of the Control/Test Bead, for air bubbles left in the chamber and liquid in the top of the device. If undissolved beads, air bubbles or liquid in the top of the device are noted, tap the device on the bench top and repeat visual inspection. Amplification and detection should be initiated within 15 minutes.
5. Insert the *illumigene* Test Device into the *illumipro-10* and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

| Sample ID | Reported Result | Interpretation |
|------------------|-----------------|---|
| Patient Specimen | POSITIVE | Sample contains <i>Bordetella pertussis</i> target DNA. |
| | NEGATIVE | No <i>Bordetella pertussis</i> DNA detected. |
| | INVALID | No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure. |
| Positive Control | POSITIVE | Valid positive control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly. |
| | NEGATIVE | Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. |
| | INVALID | No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure. |
| Negative Control | POSITIVE | Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services Department at 1-800-343-3858 (US) or your local distributor. |
| | NEGATIVE | Valid negative control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly. |
| | INVALID | No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure. |
| EMPTY WELL | NONE | No <i>illumigene</i> Test Device in the <i>illumipro-10</i> Well. OR The <i>illumigene</i> Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample. |

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control chamber that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat-treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 ± 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 ± 2 minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. *illumigene* Pertussis External Control Reagents are supplied separately (Catalog 279930). It is recommended that reactivity of each new lot and each new shipment of *illumigene* Pertussis be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The *illumigene* Pertussis test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate device must be used for each external control reagent.

EXPECTED VALUES

Overall incidence of *B. pertussis* in prospectively collected and tested specimens during the period of this study was approximately 10.0% (51/508).

LIMITATIONS OF THE PROCEDURE

1. This product can be used only with the *illumipro-10* instrument.
2. The *illumigene* Pertussis DNA assay is a qualitative assay and does not provide quantitative values or information about organism load.
3. This device has not been evaluated for monitoring treatment of *Bordetella pertussis* infections.
4. Respiratory infections can be caused by *Bordetella pertussis* as well as other pathogens. Positive results do not preclude coinfection with other respiratory pathogens.
5. The detection of nucleic acid is dependent upon proper specimen collection, handling, transportation, storage and preparation. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
Prevalence of *B. pertussis* infection will affect the test's predictive value.
7. Organism nucleic acids may persist *in vivo*, independent of organism viability. The *illumigene* Pertussis assay does not distinguish between viable and nonviable organisms.
8. As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen; (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.
9. The *illumigene* Pertussis assay targets the IS481 insertional element of the *Bordetella* genome. The IS481 insertional element is present in *B. pertussis*, *B. holmesii* and some strains of *B. bronchiseptica*.
10. Acetyl salicylic acid, as found in aspirin, produced invalid results when tested at concentrations above 5 mg/mL during *B. pertussis* strain BAA-589 Limit of Detection replicate testing.

SPECIFIC PERFORMANCE CHARACTERISTICS

The *illumigene* Pertussis DNA Amplification Assay was evaluated from December 2012 to July 2013 by independent clinical test sites representing geographically distinct regions throughout the United States. A total of 729 qualified nasopharyngeal (NP) swab specimens collected from patients suspected of respiratory infections due to *B. pertussis* were evaluated with the test device to establish performance characteristics. Specimens were leftover, de-identified specimens that had previously been submitted for routine *B. pertussis* testing. Specimens included in performance evaluation were prospective (never frozen) and retrospective (frozen prior to *illumigene* testing). The performance of *illumigene* Pertussis was compared to a Composite Reference Method that included two manufacturer validated, IS481-targeted real-time PCR assays followed by bi-directional sequencing of amplicon from PCR positive specimens. The calculated cycle threshold (Ct) value for each PCR assay was determined to be 35.6 at the assay limit of detection (PCR1: 3750 CFU/mL, PCR2: 4522 CFU/mL). Specimens were considered positive when results from either comparator PCR assay produced Ct values below 35.6 and the corresponding bi-directional sequencing result was positive. Specimens were considered to be negative when results from both comparator PCR assays produced Ct values above 35.6. A total of 508 (69.7%) prospective and 221 (30.3%) retrospective specimens were tested. Invalid results were obtained for 13 specimens (1.8%); two specimens remained invalid after repeat testing. Two prospective and two retrospective specimens produced indeterminate comparator PCR results. Table 1 summarizes *illumigene* performance characteristics; Table 2 describes *illumigene* performance characteristics by Site.

Table 1. *illumigene* Assay Performance

| Specimen Description | Positive Specimens | | | Negative Specimens | | | Invalid Results* |
|------------------------------------|----------------------------------|--------------|---------------------|----------------------------------|--------------|---------------------|------------------|
| | <i>illumigene</i> vs. Comparator | PPA | 95% CI | <i>illumigene</i> vs. Comparator | NPA | 95% CI | |
| Composite Method Comparator | | | | | | | |
| Prospective | 30/32 | 93.8% | 79.9 – 98.3% | 451/472 | 95.6% | 93.3 – 97.1% | 2 (13) |
| Retrospective | 21/21 | 100.0% | 84.5 – 100.0% | 192/198 | 97.0% | 93.5 – 98.6% | 0 |
| TOTAL | 51/53 | 96.2% | 87.2 – 99.0% | 643/670 | 96.0% | 94.2 – 97.2% | 2 (13) |

* 11/13 initial invalid specimens produced valid results upon repeat testing

Table 2. *illumigene* Assay Performance by Site

| Specimen Description | Positive Specimens | | | Negative Specimens | | | Invalid Results* |
|------------------------------------|----------------------------------|--------------|---------------------|----------------------------------|--------------|---------------------|------------------|
| | <i>illumigene</i> vs. Comparator | PPA | 95% CI | <i>illumigene</i> vs. Comparator | NPA | 95% CI | |
| Composite Method Comparator | | | | | | | |
| Site 1 | 41/43 | 95.3% | 84.5 – 98.7% | 449/470 | 95.5% | 93.3 – 97.1% | 1 (12) |
| Site 2 | 3/3 | 100.0% | 43.9 – 100.0% | 67/70 | 95.7% | 88.1 – 98.5% | 1 (1) |
| Site 3 | 0/0 | N/A | N/A | 8/8 | 100.0% | 67.6 – 100.0% | 0 |
| Site 4 | 7/7 | 100.0% | 64.6 – 100.0% | 119/122 | 97.5% | 93.0 – 99.2% | 0 |
| TOTAL | 51/53 | 96.2% | 87.2 – 99.0% | 643/670 | 96.0% | 94.2 – 97.2% | (13) |

* 11/13 initial invalid specimens produced valid results upon repeat testing

Clinical studies were conducted with multiple nasopharyngeal swab and sample elution buffer types. Sample buffers tested during clinical studies included 0.85% Saline (n=30 or 4.1%), Tris EDTA (n=8 or 1.1%) and Molecular Grade Water (n= 687 or 94.2%). All sample buffers were used in 0.5 mL volumes. Analytical studies were performed with 0.85% Saline, Tris EDTA, Phosphate Buffered Saline (PBS) and Molecular Grade Water. Analytical studies established equivalence between all sample elution buffer types.

Age information was known for 723 (99.2%) of the patients from whom samples were tested. Patient age ranged from 1 month to 88 years. Thirty-eight (5.2%) patients were less than 1 year of age; 13 (1.8%) were between 1 and 2 years old; 296 (40.6%) were between 2 and up to 12 years, 157 (21.5%) were between 12 and up to 21 years, 190 (26.0%) were above 21 but below 65 years, and the remaining 29 (4.0%) patients were above 65 years of age. There is no expectation that the *illumigene* Pertussis assay performs differently when evaluated specimens from patients in different age groups.

The study population included 413 (56.7%) female and 308 (42.2%) male patients. Gender was unknown for 8 (1.1%) patients included in the study. There is no expectation that the *illumigene* Pertussis assay performance characteristics are influenced by patient gender.

ANALYTICAL SENSITIVITY

The analytical sensitivity or Limit of Detect for the *illumigene* Pertussis assay was determined for *Bordetella pertussis* strain ATCC BAA-589, Tahoma I.

Limit of Detection was determined using 60 replicates *Bordetella pertussis* strain ATCC BAA-589 and a stated probability (eg. 95%, where 57/60 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below.

| <i>B. pertussis</i> Strain Description | CFU/mL | CFU/Test |
|--|--------|----------|
| ATCC BAA-589 (Tahoma I) | 3265 | 1.48 |

ASSAY REACTIVITY

The following *B. pertussis* strains were tested and produced positive reactions at 3265 CFU/mL or 1.48 CFU/Test with *illumigene* Pertussis: ATCC 12743; ATCC 8478; ATCC 8467; ATCC 9797; ATCC 53894; ATCC 10380; ATCC 12742; and A639. The following *B. pertussis* strains were tested and produced positive reactions at 3500 CFU/mL or 1.59 CFU/test with *illumigene* Pertussis: ATCC 51445 and ATCC BAA-1335.

REPRODUCIBILITY

Reproducibility studies were carried out by three of the four participating Clinical Sites. Blind-coded panels of 10 samples were supplied to participating laboratories. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as moderate positive samples (1.31×10^4 CFU/mL or 6 CFU/test), low positive samples (4.89×10^3 CFU/mL or 2 CFU/test); and high negative samples (8.2 CFU/mL or 0.004 CFU/test). The panel also included one negative sample, positive control and negative control. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of *illumigene* Pertussis and six *illumipro-10* instruments were used in this study. Positive and Negative Controls were tested each day of testing. The results are provided in the table below:

| Sample Type | Site 1 | | Site 2 | | Site 4 | | Total | |
|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|--------|
| | Percent Agreement | Percent Agreement | Percent Agreement | Percent Agreement | Percent Agreement | Percent Agreement | Percent Agreement | |
| Moderate Positive | 30/30 | 100.0% | 30/30 | 100.0% | 30/30 | 100.0% | 90/90 | 100.0% |
| Low Positive | 27/30 | 90.0% | 29/30 | 96.7% | 30/30 | 100.0% | 86/90 | 95.6% |
| High Negative | 26/30 | 86.7% | 23/30 | 76.7% | 29/30 | 96.7% | 78/90 | 86.7% |
| Negative | 10/10 | 100.0% | 9/10 | 90.0% | 10/10 | 100.0% | 29/30 | 96.7% |
| Negative Control | 10/10 | 100.0% | 10/10 | 100.0% | 10/10 | 100.0% | 30/30 | 100.0% |
| Positive Control | 10/10 | 100.0% | 10/10 | 100.0% | 10/10 | 100.0% | 30/30 | 100.0% |

CROSSREACTIVITY

Crossreactivity studies were performed with positive and negative nasal wash specimens inoculated with bacterial or fungal organisms to a final concentration of 1.0×10^5 CFU/mL or virus at a minimum of 1.0×10^5 TCID₅₀/mL. None of the following organisms reacted with *illumigene* Pertussis: *Acinetobacter baumannii*, *Acinetobacter calcoaceticus*, *Acinetobacter Iwoffii*, *Actinomyces odontolyticus*, *Arcanobacterium haemolyticum*, *Bacillus subtilis*, *Bacteroides fragilis*, *Bordetella avium*, *Bordetella hinzii*, *Bordetella parapertussis*, *Bordetella petrii*, *Bordetella trematum*, *Burkholderia cepacia*, *Candida albicans*, *Candida glabrata*, *Chlamydia pneumoniae*, *Chlamydia trachomatis*, *Citrobacter freundii*, *Clostridium difficile*, *Corynebacterium diphtheriae*, *Corynebacterium pseudodiphtheriticum*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* (ESBL), *Fusobacterium nucleatum*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae* (KPC), *Lactobacillus acidophilus*, *Lactobacillus plantarum*, *Legionella jordanis*, *Legionella longbeachae*, *Legionella micdadei*, *Legionella pneumophila*, *Listeria monocytogenes*, *Moraxella catarrhalis*, *Mycobacterium tuberculosis*, *Mycoplasma genitalium*, *Mycoplasma hominis*, *Mycoplasma pneumoniae*, *Neisseria cinerea*, *Neisseria elongata*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Nocardia asteroides*, *Peptostreptococcus anaerobius*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Serratia liquefaciens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus anginosus* (Group F), *Streptococcus bovis* (Group D), *Streptococcus canis* (Group G), *Streptococcus dysgalactiae* ssp. *dysgalactiae*, *Streptococcus dysgalactiae* ssp. *equisimilis*, *Streptococcus intermedius*, *Streptococcus mitis*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Streptococcus salivarius*, *Streptococcus suis*, *Ureaplasma urealyticum*, Adenovirus, Coronavirus, Coxsackievirus, Cytomegalovirus, Epstein Barr Virus, Herpes Simplex Virus 1, Herpes Simplex Virus 2, Human Metapneumovirus, Influenza A, Influenza B, Measles virus, Mumps virus, Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Respiratory syncytial virus A, Respiratory syncytial virus B, Rhinovirus.

Bordetella bronchiseptica Strain 4617 and *Bordetella holmesii* were tested at 1.0×10^6 CFU/mL and were found to react with the *illumigene* Pertussis assay.

TESTS FOR INTERFERING SUBSTANCES

The following chemical substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Acetaminophen (10 mg/mL), Advil® [ibuprofen (10 mg/mL)], Afrin® Decongestant Nasal Spray [Oxymetazoline hydrochloride (0.0005% w/v)], Albuterol Sulfate [salbutamol sulfate (1% w/v)], Aspirin (5mg/mL), Coricidin® HBP Cold/Flu Tablets [Acetaminophen (3.26 mg/mL), Chlorpheniramine maleate (0.02 mg/mL)], Diphenhydramine HCl (0.25 mg/mL), Erythromycin (2% w/v), Mupirocin (2% w/v), Petroleum Jelly [white petrolatum (1% w/v)], Robitussin® Cough+Chest Congestion DM Cough Syrup [dextromethorphan HBr (0.1 mg/mL), Guaifenesin (1.0 mg/mL)], Suphedrine PE [phenylephrine HCl (0.3 mg/mL)], Saline Nasal Spray [sodium chloride (0.0065% w/v)], Smokeless Tobacco (snuff) (1% w/v), Tobramycin (0.6 mg/mL), Vicks® VaporRub® [camphor (0.48% w/v), eucalyptus oil (0.12% w/v), menthol (0.26% w/v)].

Aspirin was found to interfere with *illumigene* Pertussis testing at concentrations greater than 5 mg/mL.

The following biological substances, at the specified saturated solvent/diluent concentrations, do not interfere with test results: Human DNA (200 ng/μL), Mucin [bovine submaxillary gland type I-S (1% w/v)], Whole blood (1% w/v).

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REFERENCES

- Nagamine, K., Hase, T., Notoni, T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002; 16:223-29.
- Mori, Y., Kitao, M., Tomita, N., Notoni, T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004; 59:145-47.
- Centers for Disease Control and Prevention. 2005 CDC Guidelines: Recommended Antimicrobial Agents for the Treatment and Postexposure Prophylaxis of Pertussis. 2005; 54(RR14); 1-16.
- Centers for Disease Control and Prevention. 2012 Final Pertussis Surveillance Report. 2013; 61(52).
- Klein NP, et al. Waning Protection after Fifth Dose of Acellular Pertussis Vaccine in Children. *N Engl J Med* 2012; 367: 1012-1019.
- Misegades LK, et al. Association of Childhood Pertussis with Receipt of 5 Doses of Pertussis Vaccine by Time Since Last Vaccine Dose, California, 2010. *JAMA* 2012; 308(20): 2126-2132.
- Andre P, et al. Comparison of Serological and Real-Time PCR Assays to Diagnose Bordetella pertussis Infection in 2007. *J Clin Micro* 2008; 46(5): 1672-1677.
- Centers for Disease Control and Prevention. Best Practices for Health Care Professionals on the use of Polymerase Chain Reaction (PCR) for Diagnosing Pertussis.
- US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
- CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11189

REV. 12/13

| | |
|---|--|
| <p>Meridian Bioscience, Inc. USA/Corporate Office 3471 River Hills Drive Cincinnati, Ohio 45244 Telephone: (513) 271-3700 Orders/Customer Service: (800) 543-1980 Technical Support Center: (800) 343-3858 Information Fax:(513) 272-5432 Ordering Fax:(513) 271-0124</p> | <p>Meridian Bioscience Europe s.a. / n.v. Rue de l'Industrie 7 1400 Nivelles Belgium Tel: +32 (0) 67 89 59 59 Fax: +32 (0) 67 89 59 58 e-mail: info.bn@meridianbioscience.eu</p> |
| | <p>Meridian Bioscience Europe b.v. Halderheiweg 6 5282 SN Boxtel The Netherlands Tel: +31 (0) 411 62 11 66 Fax: +31 (0) 411 62 48 41 e-mail: info.bn@meridianbioscience.eu</p> |
| <p>Meridian Bioscience Europe Via dell'Industria, 7 20020 Villa Cortese (MI) Italy Tel: +39 0331 433636 Fax: +39 0331 433616 e-mail: info@meridianbioscience.eu</p> | <p>Meridian Bioscience Europe France 34, rue de Ponthieu 75008 Paris France Tel: +33 (0) 1 42 56 04 40 Fax: +33 (0) 9 70 06 62 10 e-mail: info.fr@meridianbioscience.eu</p> |

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product.

Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de simbolos, Erläuterung der graphischen symbole)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.