NADAL[®] Legionella/S. pneumoniae Test (test cassette)

REF 572005

| | | | |
|--------------------------------------|----|------------------------------------|----|
| DE Gebrauchsanweisung | 2 | PL Sposób użycia | 21 |
| EN Instructions for use | 6 | PT Instruções de Utilização | 25 |
| FR Instructions d'utilisation | 9 | Symbols | 31 |
| ES Instrucciones de uso | 13 | Our Teams | 32 |
| IT Istruzioni per l'uso | 17 | | |

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14. Referenzen

1. Edelstein, P.H. and Meyer, R.D., 'Legionella Pneumonias' In Respiratory infection: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New-York: Raven Press, Ltd. 1994.
2. Marston, B.J., Lipman, H.B. and Breiman, R.F. Arch. Intern. Med., 154, 2417-2422 (1994).
3. Horwitz, M.A., Marston, B.J., Broone, C.V. and Freiman, R.F. Barbaree, J.M., Breiman, R.F. and Dufour, A.P., eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology (1993).
4. Kohler, R. B. Winn, W.C. Jr and Wheat, L.J. J. Clin. Microbiol., 20, 605-607 (1984).
5. Lindsay, D.S.J., Abraham, W.H., Christie, P., Johnston, F. and Edwards, G.F.S. J. Med. Microbiol., 53, 183-187 (2004).
6. Cloud, J.L., Carroll, K.C., Pixton, P., Erali, M. and Hillyard, D.R. J. Clin. Microbiol., 38, 1709-1712.
7. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae : virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
8. CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34): 1102-1106, 2010.
9. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi : 10.1038/nrmicro1871, 2008.
10. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285: 1729-1735, 2001.
11. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343: 1917-1924, 2000.
12. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

Rev. 1, 2017-03-07 OM/UJ

1. Intended Use

The NADAL® Legionella/S. pneumoniae Test is a lateral flow, rapid immunochromatographic test for the qualitative detection of *Legionella pneumophila* and *Streptococcus pneumoniae* antigens in human urine. The test is intended to be used as an aid in the diagnosis of *Legionella* infection (Legionnaires' disease), pneumonia, meningitis and otitis.

2. Introduction and Clinical Significance

Legionella are gram-negative bacteria associated with Pontiac fever and Legionnaires' disease leading to the infection of the upper respiratory tract and causing the pulmonary disease of pneumonia. In many cases high fever has been reported.(1) An estimated 10,000 to 25,000 cases involving *Legionella pneumophila* infections occur in the USA annually.(2, 3) About 5% to 30% of people with Legionnaires' disease die as a result. Extrapulmonary infections may be spread through the bloodstream or the lymph system infecting the heart, brain, kidney, liver and/or spleen. Common methods for the detection of *Legionella pneumophila* include culture, direct fluorescent antibody testing, indirect fluorescent antibody, enzyme immunoassays and polymerase chain reaction.(5, 6) *Streptococcus pneumoniae* is a gram-positive bacterium first isolated by Pasteur in 1881 from the saliva of a patient with rabies.(7) The chemical structure and antigenicity of the pneumococcal capsular polysaccharide and its association with virulence and role in human disease were explained over the period from 1915 to 1945. Bacteria are lancet-shaped anaerobic organisms which are spread by close person-to-person contact via respiratory droplets and cause serious diseases in humans: The 10 most common serotypes are estimated to account for about 62% of invasive disease worldwide.(8) *Streptococcus pneumoniae* colonizes upper respiratory tract tissues, causing severe pneumonia and mild/acute earache/otitis.(9) Pneumococci cause 13% to 19% of all cases of bacterial meningitis in the United States.(10) One-quarter of patients with pneumococcal meningitis also have pneumonia. The clinical symptoms are generally similar to other forms of purulent bacterial meningitis and include headache, lethargy, vomiting, irritability, fever, nuchal rigidity, cranial signs, seizures and coma. The case-fatality of pneumococcal meningitis is about 30% but can be as high as 80% among the elderly. Bacterial pneumonia accounts for 12-16% of invasive pneumococcal disease among children aged 2 years and younger, and *S. pneumoniae* has become the leading cause of bacterial meningitis among children younger than 5 years old in the United States. Antibiotic treatment is efficient even if more penicillin-resistant strains have been identified.(11) Several vaccines are available, with varying efficiency, depending on patient factors such as age, chronic illnesses or immunodeficiency. Nevertheless, vaccines have been demonstrated to provide protection against pneumococcal pneumonia.(12)

3. Test Principle

The NADAL® Legionella/S. pneumoniae Test cassette consists of a plastic housing containing two different internal strips for the detection of *Legionella* and *Streptococcus*. Coloured monoclonal *Legionella* and *Streptococcus* antibodies conjugated with colloidal gold are placed in the respective conjugate pad of each membrane. A few drops of urine

sample should be dispensed into each sample well (⇔) of the test cassette. As the test sample flows along the membranes, the labelled antibody-dye conjugate binds to the *Legionella* and *Streptococcus* antigens (when present in the sample) forming antibody-antigen complexes. These complexes bind to the polyclonal *Legionella* or *Streptococcus* antibodies in the test line region (T) of each test cassette, producing a rose-pink line. In the absence of *Legionella* or *Streptococcus*, no lines will appear in the test line regions (T). The reaction mixture continues migrating along the membranes, through the control line regions (C). Further antibody conjugates bind to the reagents in the control line regions (C), producing rose-pink lines. The control line should always appear in each control line region (C) regardless of the presence or absence of the test line, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 NADAL® Legionella/S. pneumoniae test cassettes (incl. disposable pipettes)
- 1 package insert

5. Additional Materials Required

- Timer
- Specimen collection container

6. Storage & Stability

The test should be stored in the sealed foil pouch at 4-30°C. Do not freeze the test kit. The test is stable until the expiry date stated on the package.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Read the package insert carefully before performing the test.
- Do not use test kits beyond the expiration date.
- Use only urine samples.
- Do not reuse tests.
- Do not use the test if the foil pouch is damaged.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Handle all specimens as if they contained infectious agents.
- Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Used testing materials should be discarded according to local regulations.
- Humidity and temperature can adversely affect test results.


8. Specimen Collection and Preparation

Collect a urine specimen in a clean container which is free of detergent. For optimal results, a first morning urine specimen is preferred, since this sample could contain the highest concentration of bacteria. If testing is not immediate, the specimen should be refrigerated (2-8°C) for up to 24 hours. In such cases, bring the specimen to room temperature prior to testing. If testing is delayed for more than 24 hours, the specimen should be frozen. A frozen specimen must be

completely thawed, thoroughly mixed and brought to room temperature prior to testing. Avoid repeated freezing and thawing of specimens.

9. Test Procedure

1. Bring all test components and samples to room temperature.
2. Remove the test cassette from the foil pouch.
3. Label the test cassette with the patient's name or control number.
4. Fill a pipette with urine. Holding it vertically, add 3 drops of the urine to each sample well (⇨).
5. Read test results after 10 minutes. Do not interpret test results after more than 15 minutes.

10. Result Interpretation

Positive for *Legionella pneumophila*:

One coloured line develops in each control line region (C) and one coloured line develops in the test line region (T) for *Legionella pneumophila* (LEGIO).

Positive for *Streptococcus pneumoniae*:

One coloured line develops in each control line region (C) and one coloured line develops in the test line region (T) for *Streptococcus pneumoniae* (S.PNEU).

Positive for *Legionella pneumophila* and *Streptococcus pneumoniae*:

In addition to the control lines (C), a red/pink line develops in each test line region for *Legionella pneumophila* and *Streptococcus pneumoniae*.

Negative:

One coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

Invalid:

The control lines fail to appear. Results from any test which has not produced control lines at the specified reading time must be discarded.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

12. Limitations

- The NADAL® Legionella/S. pneumoniae Test is designed to detect *Legionella* and *Streptococcus* antigens in human urine samples. A negative antigen result does not rule out the possibility of infection with *Legionella pneumophila*. The antigen concentration may fall below the detection level of the test. It is recommended that culture be used for patients with suspected pneumonia to detect the presence of causative agents other than *Legionella pneumophila*.
- If a false positive result is suspected (i.e. positive result without any clinical signs of the disease), the test must be repeated on the same urine sample after the urine sample being heated for 5 minutes at 100°C to denature urine proteins without affecting *Legionella* antigens.
- False positive results can be obtained for patients having been treated with injection of anti-lymphocyte rabbit serum. Anti-rabbit antibodies can be eliminated by heating a urine specimen without affecting *Legionella* antigens (as described in the previous point).
- The test provides a presumptive diagnosis of a *Legionella pneumophila* infection. However, this finding should not be based on clinical and radiological evidence alone. In order to make an accurate diagnosis, culture results, serology testing and antigen detection methods should be used in conjunction with clinical findings.
- Sensitive immunoassays may demonstrate positive results with specimens containing heterophilic antibodies. If the qualitative interpretation is inconsistent with clinical findings, further testing with an alternative method should be performed.
- Excretion of *Legionella* antigens may be seen as early as three days after the onset of symptoms. However, antigens may continue to be excreted for up to one year afterwards. A positive result with the NADAL® Legionella/S. pneumoniae Test may occur during a current and past infection.
- The NADAL® Legionella/S. pneumoniae Test will not differentiate between colonised and infected individuals.
- The NADAL® Legionella/S. pneumoniae Test cannot distinguish between different pneumococcal infections potentially induced by *Streptococcus pneumoniae*. Therefore, in order to make an accurate diagnosis, culture results, serology testing and antigen detection methods should be used in conjunction with clinical findings.
- As with any other diagnostic procedures, the physician should evaluate data obtained with this test in the light of other clinical information.

13. Performance Characteristics

1. Legionella

a) Sensitivity and specificity

A comparative study between the BinaxNOW® Legionella rapid test and the NADAL® Legionella/S. pneumoniae Test has been performed using a panel of 48 human urine samples.

The results are presented in the following table:

| | | BinaxNOW® Legionella rapid test (reference method) | | |
|------------------------------------------------|----------|-------------------------------------------------------|----------|-------|
| | | Positive | Negative | Total |
| NADAL® Legionella/ S. pneumoniae Test | Positive | 8 | 0 | 8 |
| | Negative | 0 | 40 | 40 |
| | Total | 8 | 40 | 48 |

From the above table, the diagnostic sensitivity of the NADAL® Legionella/S. pneumoniae Test is 100% (8/8) and its specificity is 100% (40/40) when compared to the reference method. Overall agreement is 100% (48/48).

b) Cross reactivity

The *Legionella* strains listed below were cultured and prepared to a working concentration of 1×10^8 CFU/mL in saline. Each test strain was run on the NADAL® Legionella/S. pneumoniae Test. None of the organisms listed below cross-reacted with the NADAL® Legionella/S. pneumoniae Test.

| | |
|-------------------------------------------------|----------------------------|
| <i>Legionella pneumophila</i> serogroup 2-14 | <i>Legionella anisa</i> |
| <i>Legionella bozemanii</i> | <i>Legionella dumoffii</i> |
| <i>Legionella feeleii</i> | <i>Legionella gormanii</i> |
| <i>Legionella hackeliae</i> | <i>Legionella micdadei</i> |
| <i>Legionella longbeachae</i> | |

Bacteria strains other than those of the genus *Legionella* were prepared at a working concentration of 3×10^8 CFU/mL in saline and run on the NADAL® Legionella/S. pneumoniae Test. None of the organisms listed below cross-reacted with the NADAL® Legionella/S. pneumoniae Test.

| | |
|---------------------------------|-----------------------------------|
| <i>Bacillus subtilis</i> | <i>Enterococcus faecalis</i> |
| <i>Enterococcus faecium</i> | <i>Micrococcus luteus</i> |
| <i>Staphylococcus aureus</i> | <i>Staphylococcus epidermidis</i> |
| <i>Streptococcus faecalis</i> | <i>Streptococcus mutans</i> |
| <i>Streptococcus pneumoniae</i> | <i>Streptococcus pyogenes</i> |
| <i>Streptococcus sanguis</i> | <i>Acinetobacter baumannii</i> |
| <i>Branhamella catarrhalis</i> | <i>Enterobacter cloacae</i> |
| <i>Escherichia coli</i> | <i>Haemophilus influenzae</i> |
| <i>Klebsiella oxytoca</i> | <i>Klebsiella pneumoniae</i> |
| <i>Proteus mirabilis</i> | <i>Proteus valgaris</i> |
| <i>Pseudomonas seruginosa</i> | <i>Salmonella enteritidis</i> |
| <i>Salmonella paratyphi A</i> | <i>Salmonella typhi</i> |
| <i>Serratia marcescens</i> | <i>Candida albicans</i> |
| <i>Candida glabrata</i> | <i>Candida tropicalis</i> |

c) Intra-assay reproducibility

The intra-assay reproducibility was evaluated by running 10 replicates of 4 sample dilutions with differing *Legionella* concentration. All the series showed converging and expected results. Intra-lot reproducibility performance of the NADAL® Legionella/S. pneumoniae Test is 100%.

2. Streptococcus

a) Sensitivity and specificity

A comparative study between BinaxNOW® Streptococcus pneumoniae rapid test and the NADAL® Legionella/S. pneumoniae Test has been performed using a panel of 16 fresh urine samples and 42 stored urine samples.

The results are presented in the following table:

| | | BinaxNOW® Streptococcus Pneumoniae rapid test (reference method) | | |
|------------------------------------------------|----------|------------------------------------------------------------------------|----------|-------|
| | | Positive | Negative | Total |
| NADAL® Legionella/ S. pneumoniae Test | Positive | 26 | 4 | 30 |
| | Negative | 7 | 21 | 28 |
| | Total | 33 | 25 | 58 |

From the above table, the diagnostic sensitivity of the NADAL® Legionella/S. pneumoniae Test is 79% (26/33) and its specificity is 84% (21/25) when compared to the reference method. Overall agreement is 81% (47/58) and is satisfactory due to the fact that long storage of the 42 samples (33 positive and 9 negative) could have caused degradation of *Streptococcus pneumoniae* antigens.

b) Cross-reactivity

Different strains of bacteria were tested using the NADAL® Legionella/S. pneumoniae Test in order to determine whether they might cause a non-specific reaction (cross-reactions).

The results are presented in the following table:

| Group | Result |
|---------------------------------|----------|
| Strep-Pneumo capsular types (8) | Positive |
| Strep-A, B, C, D, E, F, G & H | Negative |
| <i>Staphylococcus Aureus</i> | Negative |
| <i>Bacillus subtilis</i> | Negative |

c) Intra assay-reproducibility

The intra-assay reproducibility was evaluated by running 10 replicates of 5 sample dilutions having different *Streptococcus Pneumoniae* concentration. All the series showed converging and expected results. Intra-lot reproducibility performance of the NADAL® Legionella/S. pneumoniae Test is 100%.

14. References

- Edelstein, P.H. and Meyer, R.D., 'Legionella Pneumonias' In Respiratory infection: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New-York: Raven Press, Ltd. 1994.
- Marston, B.J., Lipman, H.B. and Breiman, R.F. Arch. Intern. Med., 154, 2417-2422 (1994).
- Horwitz, M.A., Marston, B.J., Broone, C.V. and Freiman, R.F. Barbaree, J.M., Breiman, R.F. and Dufour, A.P., eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology (1993).
- Kohler, R. B. Winn, W.C. Jr and Wheat, L.J. J. Clin. Microbiol., 20, 605-607 (1984).
- Lindsay, D.S.J., Abraham, W.H., Christie, P., Johnston, F. and Edwards, G.F.S. J. Med. Microbiol., 53, 183-187 (2004).
- Cloud, J.L., Carroll, K.C., Pixton, P., Erali, M. and Hillyard, D.R. J. Clin. Microbiol., 38, 1709-1712.
- AlonsoDeVelasco, E., Verheul, A.F., Verhoen, J., Snippe, H. Streptococcus pneumoniae : virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
- CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34): 1102-1106, 2010.
- Kadioglu, A., Welsler, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi : 10.1038/nrmicro1871, 2008.
- Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285: 1729-1735, 2001.
- Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343: 1917-1924, 2000.
- Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

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14. Bibliographie

1. Edelstein, P.H. and Meyer, R.D., 'Legionella Pneumonias' in Respiratory infection: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New-York: Raven Press, Ltd. 1994.
2. Marston, B.J., Lipman, H.B. and Breiman, R.F. Arch. Intern. Med., 154, 2417-2422 (1994).
3. Horwitz, M.A., Marston, B.J., Broone, C.V. and Freiman, R.F. Barbaree, J.M., Breiman, R.F. and Dufour, A.P., eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology (1993).
4. Kohler, R. B. Winn, W.C. Jr and Wheat, L.J. J. Clin. Microbiol., 20, 605-607 (1984).
5. Lindsay, D.S.J., Abraham, W.H., Christie, P., Johnston, F. and Edwards, G.F.S. J. Med. Microbiol., 53, 183-187 (2004).
6. Cloud, J.L., Carroll, K.C., Pixton, P., Erali, M. and Hillyard, D.R. J. Clin. Microbiol., 38, 1709-1712 .
7. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae : virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
8. CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34): 1102-1106, 2010.
9. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi : 10.1038/nrmicro1871, 2008.
10. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285: 1729-1735, 2001.
11. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343: 1917-1924, 2000.
12. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

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14. Bibliografia

1. Edelstein, P.H. and Meyer, R.D., 'Legionella Pneumonias' In Respiratory infection: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New-York: Raven Press, Ltd. 1994.
2. Marston, B.J., Lipman, H.B. and Breiman, R.F. Arch. Intern. Med., 154, 2417-2422 (1994).
3. Horwitz, M.A., Marston, B.J., Broone, C.V. and Freiman, R.F. Barbaree, J.M., Breiman, R.F. and Dufour, A.P., eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology (1993).
4. Kohler, R. B. Winn, W.C. Jr and Wheat, L.J. J. Clin. Microbiol., 20, 605-607 (1984).
5. Lindsay, D.S.J., Abraham, W.H., Christie, P., Johnston, F. and Edwards, G.F.S. J. Med. Microbiol., 53, 183-187 (2004).
6. Cloud, J.L., Carroll, K.C., Pixton, P., Erali, M. and Hillyard, D.R. J. Clin. Microbiol., 38, 1709-1712.
7. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J., Snippe, H. Streptococcus pneumoniae : virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
8. CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34): 1102-1106, 2010.
9. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi: 10.1038/nrmicro1871, 2008.
10. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285: 1729-1735, 2001.
11. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343: 1917-1924, 2000.
12. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

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
14. Referências

1. Edelstein, P.H. and Meyer, R.D., 'Legionella Pneumonias' In Respiratory infection: Diagnosis and Management, 3rd Edition, James E. Pennington, editor. New-York: Raven Press, Ltd. 1994.
2. Marston, B.J., Lipman, H.B. and Breiman, R.F. Arch. Intern. Med., 154, 2417-2422 (1994).
3. Horwitz, M.A., Marston, B.J., Broone, C.V. and Freiman, R.F. Barbaree, J.M., Breiman, R.F. and Dufour, A.P., eds. Legionella: Current status and emerging perspectives. Washington, DC American Society for Microbiology (1993).
4. Kohler, R. B. Winn, W.C. Jr and Wheat, L.J. J. Clin. Microbiol., 20, 605-607 (1984).
5. Lindsay, D.S.J., Abraham, W.H., Christie, P., Johnston, F. and Edwards, G.F.S. J. Med. Microbiol., 53, 183-187 (2004).
6. Cloud, J.L., Carroll, K.C., Pixton, P., Erali, M. and Hillyard, D.R. J. Clin. Microbiol., 38, 1709-1712 .
7. AlonsoDeVelasco, E., Verheul,A.F.,Verhoen, J, Snippe, H. Streptococcus pneumoniae : virulence factors, pathogenesis, and vaccines. Microbiological Reviews, 59 (4), 591 – 603, 1995.
8. CDC. Updated recommendations for prevention of invasive pneumococcal disease among adults using the 23-valent pneumococcal polysaccharide vaccine (PPSV23). MMWR, 59 (34): 1102-1106, 2010.
9. Kadioglu, A., Weiser, J.N., Paton, J.C., and Andrew, P.W. The role of Streptococcus pneumoniae virulence factors in host respiratory colonization and disease. Nature Reviews, Microbiology, 6 (4) doi : 10.1038/nrmicro1871, 2008.
10. Robinson, K.A., Baughman, W., Rothrock, G. Epidemiology of invasive Streptococcus pneumoniae infections in the United States, 1995-1998. Opportunities for prevention in the conjugate vaccine era. JAMA, 285: 1729-1735, 2001.
11. Whitney, C.G., Farley, M.M., Hadler, J. et al. Increasing prevalence of multidrug-resistant Streptococcus pneumoniae in the United States. N. Engl. J. Med., 343: 1917-1924, 2000.
12. Tsai, C.J., Griffin, M.R., Nuorti, J.P. et al. Changing epidemiology of pneumococcal meningitis after the introduction of pneumococcal conjugate vaccine in the United States. Clin. Infect. Dis., 46 : 1664-1672, 2008.

Rev. 1, 2017-03-07 AL

| Symbol | Deutsch | English | Français | Español | Italiano | Polski |
|----------------------------------------------------------------------------------|------------------------------|-------------------------------------------|--------------------------------------------------|-----------------------------------------------------|------------------------------------------------|--------------------------------------|
|  | CE Konformitätszeichen | CE marking of conformity | Conforme aux normes européennes | Conformidad europea | Conformità europea | Znak zgodności CE |
|  | Gebrauchsanweisung beachten | Consult instructions for use | Consulter la notice d'utilisation | Consúltense las instrucciones de uso | Consultare le istruzioni per l'uso | Przestrzegać instrukcji obsługi |
|  | <i>in-vitro</i> -Diagnostika | <i>in-vitro</i> diagnostic medical device | Dispositif médical de diagnostic <i>in-vitro</i> | Producto sanitario para diagnóstico <i>in-vitro</i> | Dispositivo medico-diagnostico <i>in-vitro</i> | Tylko do diagnostyki <i>in-vitro</i> |
|  | Temperaturbegrenzung | Temperature limitation | Limites de température | Limite de temperatura | Limiti di temperatura | Temperatura przechowywania |
|  | Chargenbezeichnung | Batch code | Numéro de lot | Código de lote | Codice lotto | Numer serii |
|  | Nicht zur Wiederverwendung | Do not reuse | Ne pas réutiliser | No reutilizar | Non riutilizzare | Tylko do jednorazowego użytku |
|  | Verwendbar bis | Use by | Utiliser jusqu'au | Fecha de caducidad | Utilizzare entro | Data ważności |
|  | Bestellnummer | Catalogue Number | Référence du catalogue | Número de catálogo | Riferimento di Catalogo | Numer katalogowy |
|  | Hersteller | Manufacturer | Fabricant | Fabricante | Fabbricante | Producent |
|  | Ausreichend für <n> Ansätze | Sufficient for <n> tests | Suffisant pour "n" tests | Suficiente para <n> utilizaciones | Sufficiente per "n" saggi | Wystarczający na <n> Powtórzeń |

| Symbol | Português | Český | Suomi | Svenskt | Nederlands | Dansk | Norsk |
|------------------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------|--------------------------------------------------------------------|---------------------------------------------------------------|-----------------------------------------------------|--------------------------------------------------|--------------------------------------------|
|  | Conformidade com as normas europeias | CE certifikát | CE-merkitty | CE-märkning | CE-markering | CE-mærkning | CE standardisert |
|  | Consultar as instruções de utilização | Viz návod k použití | Katso käyttöohjetta | Läs bruksanvisningen | Raadpleeg de gebruiksaanwijzing | Se brugsanvisningen | Les bruksanvisning nøye |
|  | Dispositivo médico para diagnóstico <i>in-vitro</i> | Diagnostický zdravotnický prostředek <i>in-vitro</i> | <i>in-vitro</i> - diagnostiikkaan tarkoitettu lääkinnällinen laite | Medicinteknisk produkt avsedd för <i>in-vitro</i> -diagnostik | Medisch hulpmiddel voor <i>in-vitro</i> diagnostiek | Medicinsk udstyr til <i>in-vitro</i> -diagnostik | <i>in-vitro</i> diagnostic medisinsk enhet |
|  | Limites de temperatura | Teplotní omezení | Lämpötilarajat | Temperaturbegränsning | Temperatuurlimiet | Temperaturbegrænsning | Temperaturbegrensning |
|  | Código do lote | Kód šarže | Eräkoodi | Satsnummer | Code van de partij | Batchkode | Merking |
|  | Não reutilizar | Pro jednorázové použití | Kertakäyttöinen | Får inte återanvändas | Niet opnieuw gebruiken | Må ikke genbruges | Må ikke brukes om igjen |
|  | Prazo de validade | Spotřebuje do | Käytettävä viimeistään | Används före | Houdbaar tot | Udløbsdato | Tidtaking |
|  | Número de catálogo | Katalogov číslo | Luettelonumero | Listnummer | Catalogus nummer | Best il l ingsnummer | Katalog nummer |
|  | Fabricante | Výrobce | Valmistaja | Tillverkare | Fabrikant | Fabrikant | Produsent |
|  | Suficiente para <n> test | Dostačuje pro <n> testů | Lukumäärä <n> test | Räcker till <n> test | Voldoende voor <n> test | Tilstrækkeligt til <n> test | Tilstrækkelig for <n> tester |

Our Teams

Germany:

Regensburg

Tel: +49 941 290 10-0
Fax: +49 941 290 10-50

Moers

Tel: +49 2841 99820-0
Fax: +49 2841 99820-1

Austria:

Tel: +49 941 290 10-29
Free Tel: 0800 291 565
Fax: +49 290 10-50
Free Fax: 0800 298 197

UK & Ireland:

Tel: +49 941 290 10-18
Free Tel – UK: 0808 234 1237
Free Tel – IRE: 1800 555 080
Fax: +49 290 10-50

France:

France Tel: 0800 915 240
France Fax: 0800 909 493

Switzerland

Swiss Tel: 0800 564 720
Swiss Fax: 0800 837 476

Belgium

Belgium Tel: 0800 718 82
Belgium Fax: 0800 747 07

Luxembourg

Lux. Tel: 800 211 16
Lux. Fax: 800 261 79

Spain:

Tel: +49 941 290 10-759
Free Tel: 900 938 315
Fax: +49 941 290 10-50
Free Fax: 900 984 992

Italy:

Tel: +49 941 290 10-34
Fax: +49 941 290 10-50

Poland:

Tel: +49 941 290 10-44
Free Tel: 00 800 491 15 95
Fax: +49 941 290 10-50
Free Fax: 00 800 491 15 94

Portugal:

Tel: +49 941 290 10-735
Tel. Verde: 800 849 230
Fax: +49 941 290 10-50
Fax Verde: 800 849 229

Netherlands:

Tel: +31 30 75 600
Free Tel: 0800 0222 890
Fax: +31 70 30 30 775
Free Fax: 0800 024 9519

Nordic countries (Finland, Norway, Sweden, Denmark):

Tel: +31 703075 607
Free Tel: +45 80 88 87 53
Tax: +31 703030 775

Laboratory Diagnostics Team:

Tel: +49 941 290 10-40
Fax: +49 941 290 10-50

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Carl-Zeiss-Strasse 12 • 47445 Moers • Germany
www.nal-vonminden.com • info@nal-vonminden.com
Fon: +49 2841 99820-0 • Fax: +49 2841 99820-1

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nal von minden GmbH