**ImmunoCard**  
**STAT**  **NOROVIRUS**

REF 757020

IVD

ENGLISH p. 1

RAPID IMMUNOASSAY FOR THE DETECTION OF NOROVIRUS ANTIGENS  
IN HUMAN STOOL

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# ImmunoCard STAT! NOROVIRUS

**REF** Catalogue Number 757020

**IVD** *In Vitro* Diagnostic Medical Device

## RAPID IMMUNOASSAY FOR THE DETECTION OF NOROVIRUS ANTIGENS IN HUMAN STOOL

### INTENDED USE

The ImmunoCard STAT! Norovirus Immunoassay is a rapid *in vitro* qualitative test for the detection of Norovirus antigens in human stool. The test results are intended to aid in the diagnosis of Norovirus associated gastroenteritis.

### SUMMARY AND EXPLANATION OF THE TEST

Noroviruses are a major cause of acute gastroenteritis, they are involved in outbreaks in closed and crowded environments, such as schools, hospitals, nursing homes and cruise ships<sup>1,2,3</sup>. In children Noroviruses are the second most common cause of gastroenteritis after Rotavirus<sup>4</sup>. Typically, people with Norovirus infection develop diarrhoea and abdominal pain and begin to vomit within 24 to 48 hours of exposure. Norovirus symptoms last one or two days, and most people recover completely without treatment. However, in some people especially infants, older adults and people with underlying disease, vomiting and diarrhoea can be severely dehydrating and require medical attention. Noroviruses are highly contagious, commonly spread through food or water contaminated by faecal matter. The infection can also be acquired through close contact with an infected person. Gastroenteritis outbreaks caused by Norovirus are more frequent than outbreaks due to bacterial pathogens, can be difficult to control and long-lasting because Norovirus is easily transmitted due to the low infectious dose and persistent in the environment resisting disinfection and chlorination. Based on their capsid gene (VP1) sequence, human Norovirus strains can be grouped into two genogroups (GI and GII)<sup>5</sup>.

### BIOLOGICAL PRINCIPLES

The ImmunoCard STAT! Norovirus assay detects the presence of Norovirus antigens in stool. Patient specimen collected with the provided swab is diluted in the Sample Diluent Tube. The suspension is mixed and 3 drops are added to the sample port of the device. The sample mobilizes latex particles (blue) coated with anti-NV-GI (Norovirus Genogroup I) and anti-NV-GII (Norovirus Genogroup II) monoclonal antibodies and migrates along the membrane through the **Test** [immobilized anti-NV-GI, anti-NV-GII and anti-NV (reacting with both GI and GII) monoclonal antibodies] and **Control** (immobilized rabbit anti-mouse immunoglobulin) lines. After 15 minutes, the **Test** and **Control** zones are observed for the presence of blue lines across the membrane surface. If Norovirus is present in the sample, it forms a complex with the capture antibody and the monoclonal antibodies-latex conjugate which can be seen visually as a blue line in the **Test** zone. No blue line in the **Test** zone indicates a negative result. The **Control** line serves as a procedural control to assure that the sample has migrated the appropriate distance along the membrane and the test has worked properly.

### REAGENTS/MATERIALS PROVIDED

1. **Test Cards:** individually foil pouched devices containing immobilized anti-NV-GI (Genogroup I), anti-NV-GII (Genogroup II) and anti-NV (Genogroup I and II) monoclonal antibodies (**Test** zone), immobilized rabbit anti-mouse immunoglobulin (**Control** zone) and embedded anti-NV-GI and anti-NV-GII monoclonal antibodies conjugated to latex particles.
2. **Sample Diluent/Negative Control Tubes:** buffer solution containing surfactant and 0.08% sodium azide as a preservative (400 µl/tube).
3. **Swabs:** cotton/rayon unsterilized applicator sticks for stool collection.
4. **Dropper Caps with Filter:** plastic caps with dropper and integrated filter for sample filtration and dispensing  
The maximum number of tests obtained from this kit is listed on the outer box.

### MATERIALS NOT PROVIDED

1. **Positive Control** (Meridian Catalogue # 757001, 1 swab): cotton/rayon swab with recombinant non infectious Norovirus antigens.
2. **Timer**

## PRECAUTIONS


1. All reagents are for *in vitro* diagnostic use only and are ready to use.
2. Reagent concentration, incubation times and temperatures (15-30°C) have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications. Once the assay has been started, complete all subsequent steps without interruption.
3. Stool samples must be mixed thoroughly to ensure a representative sample prior to sampling. Do not use stools that have dried out.
4. For faecal sample collection be sure to use the cotton/rayon swabs provided within the ImmunoCard STAT Norovirus kit. These swabs have been specifically designed for stool collection and use with this kit. Do not use these swabs for any other application. Use the swab to collect stool excreted from the body. Do not use rectal swabs with this test.
5. Dilution of stool as described is important. Over-inoculation of stool into the Sample Diluent may cause clogging of the filter in the Dropper Cap and not allow the sample application to the Test Card. Excessive amount of collected sample may restrict movement within the Test Card so as to produce an invalid result. If clogging of the filter in the Dropper Cap or clogging of the membrane in the Test Card occurs, start over from specimen collection in a new Sample Diluent Tube and be sure to collect the indicated amount of specimen.
6. After adding the stool to the Sample Diluent Tube attach the Dropper cap firmly and test immediately without delay.
7. Dark stools may colour the Test Card membrane and make results difficult to read.
8. Patient specimens and used ImmunoCard STAT! Norovirus Test Cards may contain infectious agents and should be handled and disposed of as potentially biohazardous.
9. Store kits components at 2-30°C, avoid exposure to direct sunlight. Do not freeze.
10. Do not use the kit after its expiration date. Do not interchange reagents from different kit lot numbers.
11. Allow kit components and specimens to reach the room temperature (15-30°C) before performing a test, as cold reagents and/or specimens may decrease assay sensitivity. Reagents may take 20-30 minutes to warm following refrigeration.
12. Before performing the test verify that there is no damage or abnormality with the Swabs, Test Cards, Sample Diluent Tubes, Dropper Caps with Filter. Do not use if they are damaged.
13. If the Sample Diluent shows air bubbles at the bottom of the tube, tap the tube on a bench to eliminate them before removing the aluminium sticker.
14. When applying the specimen to the Test Card hold the Sample Diluent Tube vertically to insure proper drop size delivery.
15. Test instructions should be thoroughly read before performing any testing.

## SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. The kit components must be used before the expiration date. Store the kit at 2-30°C.

## PROCEDURAL NOTES

1. The ImmunoCard STAT! Norovirus Test Card is shown in the diagram below.
2. Batch processing of samples is possible provided correct incubation time is maintained for each device.
3. The **Control** zone of each device is a procedural control to test for proper specimen and reagent flow characteristics as well as reagent performance.

## REAGENT PREPARATION

1. All reagents are ready to use (no further dilution is required).
2. When the kit is stored in a refrigerator, allow kit components to reach room temperature (15-30°C) prior to use.
3. After removing the Test Card from its pouch or removing the aluminium sticker from the Sample Diluent Tube use them immediately.
4. Depending on the number of specimens to be tested, prepare the appropriate number of Swabs needed to collect the specimens, Sample Diluent Tubes, Dropper Caps with Filter and Test Cards.

## SPECIMEN COLLECTION AND HANDLING


For the best results, specimens should be collected as soon as possible after onset of symptoms, Norovirus concentration in stool peaks at 1-4 days after the onset of the symptoms. DO NOT USE stool in transport media or preservatives. Stool specimens may be collected by the method routinely utilized by the laboratory, providing no dilution of the specimen occurs. The sample may be collected into a clean, empty container. The specimen should be tested as soon as possible after collection, but may be stored up to 24 hours at 2-8°C prior to testing. If testing cannot be performed within this time frame, specimen should be frozen in a non-defrosting freezer (-20°C or lower) until tested.

## TEST PROCEDURE

*This test should be performed by qualified personnel per applicable local regulatory requirements.*

**Note: All components must be at 15-30°C prior to use.**

1. Remove the aluminium seal from the Sample Diluent Tube.
2. Sample Dilution: collect the stool specimen with the Swab supplied with the kit as illustrated in the diagram below. Swabs from other suppliers have not been validated.
  - **Watery stool:** collect thoroughly soaking the swab to cover the entire cotton/rayon surface.
  - **Solid or Loose stool:** collect an amount to lightly cover the cotton/rayon swab surface and rotate the swab few times.
  - The stool should be absorbed by the cotton/rayon swab as in the diagram.

3. Immediately place the swab into the Sample Diluent Tube.
4. Vigorously mix the solution by rotating the swab against the side of the Sample Diluent Tube.
5. Release as much liquid as possible from the swab by squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
6. Firmly attach the Dropper Cap with Filter to the Sample Diluent Tube.
7. Remove the Test Card from the pouch. Label appropriately. Use one card per sample or control.
8. Turn the tube upside down and dispense 3 drops of diluted specimen to **Sample** port (indicated by an arrow) of the Test Card.
9. Incubate 15 minutes at 15-30°C.  
**NOTE:** During the 15 minutes incubation, diluted specimen must move past the **Control** zone.
10. At the end of the incubation period immediately read **Control** and **Test** zones for the presence or absence of a blue line.

## EXTERNAL CONTROL TESTS

**Positive Control Swab** (Meridian Catalogue # 757001) and **Sample Diluent/Negative Control** are intended for external quality control testing as reported in the section QUALITY CONTROL.

**Positive Control Test:** after removing the Positive Control Swab from the pouch immediately place it into a Sample Diluent Tube. Continue as you would for a patient sample, as instructed at point 4 of the TEST PROCEDURE section.

**Negative Control Test:** remove the aluminium seal from a Sample Diluent Tube. Continue as you would for a patient sample, as instructed at point 6 of the TEST PROCEDURE section.

The Positive Control should yield visually detectable blue **Test** and blue **Control** lines.

The Negative Control should yield a visually detectable blue **Control** line. No **Test** line should be present.

## INTERPRETATION OF RESULTS

**Positive Test Result:** visually detectable blue **Test** and **Control** lines. The color of the Test line may vary from weak to strong blue. A positive result indicates the presence of Norovirus antigens.

**Negative Test Result:** visually detectable blue **Control** line. No **Test** line present. A negative result indicates that Norovirus antigens are absent or below the level of detection.

**Invalid Test Result:** no visually detectable blue **Control** line, with or without a visually detectable blue **Test** line.

Invalid test results may be due to a Reagent/Card problem, a procedural error, or over-inoculation of stool into Sample Diluent during specimen dilution. Re-dilute stool and repeat test.

## QUALITY CONTROL

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing or leakage. Do not use contaminated or deteriorated reagents.

### 1. Internal Control

An internal procedural control is included within each Test Card and is therefore evaluated with each test. A blue Control line serves as an internal procedural control and indicates the test has been performed correctly, that proper flow occurred and the test reagents were reactive at the time of use.

### 2. External Controls

External controls (Positive and Negative) are intended to monitor reagents reactivity and should be tested according to the requirements of the laboratory or those of applicable local, state or national accrediting agencies. To this purpose the use of Positive Control (Meridian Catalogue # 757001) and Negative Control (Sample Diluent Tube) is recommended. The Positive and Negative Controls are best tested upon receipt of the kit. Thereafter, the procedural control, performed with each test, is sufficient to verify reagent integrity and assay reactivity.

Failure of the controls to produce the expected results can mean that one of the reagents or components is no longer reactive at the time of use, the test was not performed correctly or that reagents or samples were not added. See section EXTERNAL CONTROL TESTS for instructions on performing these control tests. Positive and negative reactions are described in the section INTERPRETATION OF RESULTS. Negative Control reagent should produce a blue Control Line without any Test Lines. Positive Control contains Norovirus antigens and is expected to produce a blue Control Line and a blue Test Line.

If the expected control reactions are not observed, repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact your local distributor for technical assistance.

## EXPECTED VALUES

The ImmunoCard STAT! Norovirus test detects the presence of Norovirus antigens in stool. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on patient age, geographic location, season, method of specimen collection, handling and transportation, general health environment of the patient population under study and Norovirus genotype circulating. Noroviruses are recognized as an important cause of acute non bacterial outbreaks of gastroenteritis in all age groups. In children Noroviruses are the second most common cause of gastroenteritis after Rotavirus. Outbreaks of Norovirus have been reported to occur in a variety of settings, especially in semi closed communities such as families, schools, elderly people's homes, hospitals, hotels and cruise ships. The virus has also been identified in sporadic cases of gastroenteritis. While most of the outbreaks have a seasonal pattern, sporadic cases usually occur throughout the year<sup>6</sup>.

## LIMITATIONS OF THE PROCEDURE

1. ImmunoCard STAT! Norovirus is specifically designed to detect Norovirus antigens in stool. As with all *in vitro* diagnostic procedures, test results should be interpreted by a physician in conjunction with other clinical information.
2. A positive result does not exclude the presence of other infective organisms.
3. Do not test meconium stool and stool collected from neonates (one month or younger) since this may result in a false positive result.
4. Failure to add sufficient stool to the Sample Diluent may result in a falsely negative test result. Addition of too much stool may result in invalid test results due to the inhibition of proper sample flow.
5. Stool collected after enemas are not recommended as their performance characteristics have not been evaluated and their use could affect the test results. Stool sample obtained after eating food containing gelatinizing agents to assist swallowing could interfere with the test and are not recommended.

## SPECIFIC PERFORMANCE CHARACTERISTICS

The ImmunoCard STAT! Norovirus test was evaluated on 196 faecal specimens from healthy individuals and patients with symptoms of gastroenteritis. All specimens were tested by ImmunoCard STAT! Norovirus and compared to RT-PCR as reference test. A total of 98 stool samples resulted positive with RT-PCR for the presence of Norovirus and 98 stool samples were negative by RT-PCR. The test results are reported in the table below:

Comparison of ImmunoCard STAT! Norovirus to RT-PCR (n=196)

		RT-PCR	
		+	-
ImmunoCard STAT! Norovirus	+	80	3
	-	18	95
	<b>Total</b>	<b>98</b>	<b>98</b>

Sensitivity: 81.6% (95% confidence interval 73.9%-89.3%)

Specificity: 96.9% (95% confidence interval 93.5%-100%)

Negative Predictive Value (NPV): 84.4% (95% confidence interval 77.2%-91.6%)

Positive Predictive Value (PPV): 97.0% (95% confidence interval 93.6%-100%)

## REPRODUCIBILITY

Within-run reproducibility has been evaluated. One strong and one weak positive control for Norovirus GI and one strong and one weak positive control for Norovirus GII have been tested with ImmunoCard STAT! Norovirus, each in 5 simultaneous tests. The **Test** line of the ImmunoCard STAT! Norovirus device showed 100% reproducibility with both GI and GII strong and weak positive specimens in all tests. One Negative control specimen has been tested in 5 simultaneous tests and all test results were confirmed negative.

## CROSSREACTIVITY

The ImmunoCard STAT! Norovirus test was evaluated for specificity using the reference bacterial and viral strains listed below.

### Bacteria

Cross reactivity was not observed with the following diarrhoea-causing bacteria, intestinal flora bacteria, and other bacteria causing nosocomial or opportunistic infections ( $1 \times 10^8$  CFU/mL):

<i>Bacillus cereus</i>	<i>Escherichia coli</i> (O6)	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter coli</i>	<i>Escherichia coli</i> (O78)	<i>Salmonella enteritidis</i>
<i>Campylobacter jejuni</i>	<i>Escherichia coli</i> (O114)	<i>Salmonella typhimurium</i>
<i>Citrobacter freundii</i>	<i>Escherichia coli</i> (O126)	<i>Shigella flexneri</i>
<i>Clostridium perfringens</i>	<i>Listeria monocytogenes</i>	<i>Shigella sonnei</i>
<i>Enterococcus faecalis</i>	<i>Proteus mirabilis</i>	<i>Vibrio cholerae</i>
		<i>Vibrio parahaemolyticus</i>

### Viruses

Cross reactivity was not observed with the following diarrhoea-causing viruses:

Adenovirus type 40 ( $1 \times 10^{6.5}$ TCID <sub>50</sub> /mL)
Adenovirus type 41 ( $1 \times 10^{5.5}$ TCID <sub>50</sub> /mL)
Rotavirus A ( $1 \times 10^{5.25}$ TCID <sub>50</sub> /mL)

## TESTS FOR INTERFERING SUBSTANCES

The following substances were found to have no effect on results when present in stool at the concentrations indicated.

1. Haemoglobin: tests with added haemoglobin simulating bloody stool showed that a concentration up to 0.5 g/dL had no effect on test results.
2. Intrafat® Injection: this commercially available fat emulsion (soybean oil 10%, 2.5% glycerol, 1.2% egg phosphatide) added to stool to simulate fatty stool showed that a concentration up to 3% had no effect on test results.

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1. Rooney RM, Kramer EH, Mantha S, Nichols G, Bartram JK, Farber JM, Benembarek PK. A review of outbreaks of food borne diseases associated with passengers ships: evidence for risk management. Public Health Reports. 2004;119: 427-34.
2. Johansson PJH, Torvén M, Hammarlund A, Björne U, Hedlund K, Svensson L. Food-borne outbreak of gastroenteritis associated with Genogroup I Calicivirus. J. Clin. Microbiol. 2002; 40: 794-98.
3. Lopman BA, Reacher MH, Vipond IB, Hill D, Perry C, Halladay T, Brown DW, Edmunds WJ, Sarangi J. Epidemiology and cost of nosocomial gastroenteritis, Avon, England, 2002-2003. Emerg Infect Dis. 2004;10:1827-34.
4. Cunliffe NA, Booth JA, Elliot C, Lowe SJ, Sopwith W, Kitchin N, et al. Healthcare-associated viral gastroenteritis among children in a large pediatric Hospital, United Kingdom. Emerg Infect Dis. 2010; 16: 55-62.
5. Hansman GS, Natori K, Shirato-Horikoshi H, Ogawa S, Oka T, Katayama K, Tanaka T, Miyoshi T, Sakae K, Kobayashi S, Shinohara M, Uchida K, Sakurai N, Shinozaki K, Okada M, Seto Y, Kamata K, Nagata N, Tanaka K, Miyamura T, Takeda N. Genetic and antigenic diversity among noroviruses. J Gen Virol. 2006;87:909-19.
6. Ike Molecular Epidemiology of Norovirus in Outbreaks of Gastroenteritis in Southwest Germany from 2001 to 2004. J Clin Microbiol, 2006, 44:1262-67.

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You may see one or more of these symbols on the labeling/packaging of this product.

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**Meridian Bioscience, Inc.**