# illumigene™ *C. difficile*

## DNA Amplification Assay for the Detection of Cytotoxigenic *C. difficile* in Stool Specimens

**REF** 280050

**IVD** In vitro diagnostic medical device

### INTENDED USE

The *illumigene C. difficile* DNA amplification assay, performed on the *illumipro-10*, is a qualitative in vitro diagnostic test for the direct detection of toxigenic *C. difficile* in human stool specimens from patients suspected of having *Clostridium difficile*-associated disease (CDAD).

The *illumigene C. difficile* assay utilizes loop-mediated isothermal DNA amplification (LAMP)<sup>1,2</sup> technology to detect the pathogenicity locus (PaLoc)<sup>3</sup> of toxigenic *Clostridium difficile*. The *Clostridium difficile* PaLoc is a gene segment present in all known toxigenic *C. difficile* strains. The *C. difficile* PaLoc codes for both the Toxin A gene (*tcdA*) and the Toxin B gene (*tcdB*), has conserved border regions, and is found at the same site on the *C. difficile* genome for all toxigenic strains.<sup>3</sup> The *illumigene C. difficile* assay detects the PaLoc by targeting a partial DNA fragment on the Toxin A gene. The *tcdA* target region was selected as an intact region remaining in all known A+B+ and A-B+ toxinotypes.

*illumigene C. difficile* is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

### SUMMARY AND EXPLANATION OF THE TEST

The *illumigene C. difficile* DNA molecular assay is based on loop-mediated amplification technology, which uses specifically designed primers to the PaLoc pathogenicity locus to provide for specific and continuous isothermal DNA amplification. A by-product of this amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. This presence of turbidity signifies a positive reaction while the absence of turbidity represents a negative reaction. The *illumigene C. difficile* assay contains primers that specifically amplify a 204 bp region of the conserved 5' sequence of the *tcdA* gene within the PaLoc of toxigenic *C. difficile* in diarrheal stool samples from patients suspected of having *C. difficile*-associated disease (CDAD). The results of the assay are determined using the Meridian *illumipro-10* Incubator / Reader.

### BIOLOGICAL PRINCIPLES

Toxigenic *Clostridium difficile* is a major cause of antibiotic associated diarrhea and colitis and is the causative agent for virtually all cases of pseudomembranous colitis. Although about 2% of normal healthy adults are colonized with *C. difficile*, many patients acquire this organism through nosocomial infection. Exposure to most antibiotics is thought to allow proliferation of toxigenic *C. difficile* by disrupting the normal intestinal flora. Two large toxin proteins (TcdA [or toxin A] and TcdB [toxin B]) are thought to be the primary virulence factors of *C. difficile*. These toxins are encoded by two separate genes, named *tcdA* and *tcdB*, respectively. Together, with three additional genes, they form a 19.6 kb pathogenicity locus called PaLoc.

### REAGENTS/MATERIALS PROVIDED

**The maximum number of tests obtained from this test kit is listed on the outer box.**

- illumigene C. difficile* Sample Preparation Apparatus:** Sampling unit consisting of sample preparation chamber, dropper tip, cap and Sample Dilution Buffer (Phosphate Buffered Saline and formalin treated *Staphylococcus aureus*, with sodium azide (0.09%) as a preservative)
- illumigene* Reaction Buffer:** Tris-buffered solution containing sodium azide (0.09%) as a preservative
- illumigene C. difficile* Test Device:** Two separate chambers containing dry reagent lysospheres comprised of DNA polymerase, Deoxyribonucleoside Triphosphate (dNTPs), and either *C. difficile* specific primers (TEST Chamber) or *S. aureus* primers (CONTROL Chamber)
- Sample Collection Brushes**
- illumigene* Extraction Tubes**

### MATERIALS PROVIDED SEPARATELY

*illumigene* External Control Kit, Catalog Number: 279920

### MATERIALS NOT PROVIDED

- Disposable latex gloves, powder free
- DNase/RNase free, aerosol resistant pipette tips
- Extraction Control (known positive clinical sample or a negative sample spiked with toxigenic *C. difficile*)

### EQUIPMENT NOT PROVIDED

- Dry-bath with 12 mm heat block capable of 95 C
- Vortex mixer
- Interval timer
- Micropipette capable of dispensing 50 µL
- Micropipette capable of dispensing 200 µL
- illumipro-10*

### PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Follow Biosafety Level 2 and Good Laboratory practices during testing.<sup>4</sup> Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.

- Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
- Quality Control Programs for Molecular Testing Laboratories should be employed.<sup>5</sup>
- illumigene* Sample Dilution Buffer contains formalin inactivated-organisms and should be treated as potentially infectious.
- The *illumigene C. difficile* Test Device contains lyophilized reagents. The protective pouch should not be open until ready to perform the assay.
- The *illumigene C. difficile* Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
- Dispose of used *illumigene* Test Devices immediately after processing, leaving the device latch securely in place. Opening the device after amplification may result in contamination of the test area with amplification product.

### SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-27 C.

### SPECIMEN COLLECTION AND PREPARATION

**Sample type:** Unformed samples indicative of CDAD.

**Human stool samples, unpreserved:** Samples should be transported and stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 24 hours at 21-27 C or 5 days at 2-8 C. Samples that will not be tested within these times should be frozen immediately upon receipt and stored at ≤ -20 C until tested. Specimens may be frozen and thawed once.

**Human stool samples, preserved in Cary-Blair-based media:** Samples should be transported and stored at 2-8 C prior to testing. Samples should be tested as soon as possible, but may be held up to 5 days at 2-8 C. Samples that will not be tested within this time should be frozen immediately upon receipt and stored at ≤ -20 C until tested. Specimens may be frozen and thawed once.

### REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-27 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

### SPECIMEN PREPARATION

**NOTE:** Ensure that the *illumipro-10* instrument is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

- Mix stool sample thoroughly.
- Collect mixed sample specimen using Sample Collection Brush.
  - Liquid Stool: Immerse Sample Collection Brush completely into specimen.
  - Semi-solid Stool: Rotate Sample Collection Brush over specimen surface, lightly coating the surfaces of the brush bristles. Approximately one-half of the brush should be coated; over-collection of stool may lead to clogging of the sample collection device.
- Add the Sample Collection Brush to the *illumigene C. difficile* Sample Preparation Apparatus containing Sample Diluent and secure the cap. For stool in Cary-Blair-based media, transfer 200 µL of specimen to the *illumigene C. difficile* Sample Preparation Apparatus and secure the cap. Vortex the Sample Preparation Apparatus for a minimum of 10 seconds. Specimen in *illumigene C. difficile* Sample Preparation Apparatus prior to extraction may be held at 2-27 C for up to 24 hours prior to testing.
- Remove the tip cap from the Sample Preparation Apparatus and squeeze five to ten drops of sample into a clean *illumigene* Extraction Tube.
- Repeat Sample Preparation Steps for all samples to be processed.
- Heat the *illumigene* Extraction Tubes containing mixed sample in a dry-bath/heat block at 95 C for 10 +/- 2 minutes.
- Remove the tubes from the dry-bath/heat block and vortex for approximately 10 seconds. Extracted Samples may be held at 2-8 C for 4 hours OR may be frozen and stored at ≤ -20 C for 1 day prior to addition to Reaction Buffer. Extracted Samples may be frozen and thawed once.

### TEST PROCEDURE

**This test should be performed per applicable local, state, or federal regulations or accrediting agencies.**

**NOTE:** A maximum of 10 samples can be processed in a single *illumipro-10* run.

- Transfer 50 µL of extracted sample to an appropriately labeled *illumigene* Reaction Buffer tube.
- Vortex the Reaction Buffer Tube containing extracted sample for approximately 10 seconds.
- Repeat steps 1 and 2 for all the samples to be analyzed before proceeding.
- Remove 1 *illumigene C. difficile* Test Device from its protective pouch per sample. Carefully open the device, holding the tubes such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
- Using a new pipette tip, transfer 50 µL from the Reaction Buffer tube containing extracted sample to the TEST chamber (White Bead) of the *illumigene* Test Device. Do not insert air bubbles. Using a new pipette tip, transfer 50 µL from the Reaction Buffer tube containing extracted sample to the CONTROL chamber (Yellow Bead) of the *illumigene* Test Device. Do not insert air bubbles. Close the *illumigene* Test Device and fasten the latch securely.
- Tap device on the bench top to mix and to remove air bubbles. Carefully examine the reaction tubes to ensure that there are no air bubbles left in the tube.
- Insert the *illumigene* Test Device into the *illumipro-10* and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

## INTERPRETATION OF RESULTS

| Sample ID          | Reported Result | Interpretation  |
|--------------------|-----------------|---|
| Patient Specimen   | POSITIVE        | Sample contains toxigenic <i>C. difficile</i> strain with the pathogen locus (PaLoc).   |
|                    | NEGATIVE        | No toxigenic <i>C. difficile</i> detected.  |
|                    | INVALID         | <b>No reportable result. Repeat the test using the original stool sample.</b><br>Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.   |
| Extraction Control | POSITIVE        | Valid Extraction Control. Extractions performed correctly.  |
|                    | NEGATIVE        | <b>Extraction Failure. Do not report test result. Repeat the extraction process for all stool samples.</b>  |
|                    | INVALID         | <b>No reportable result. Do not report test result. Repeat the extraction process for all stool samples.</b><br>Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.                            |
| Positive Control   | POSITIVE        | Valid positive control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.  |
|                    | NEGATIVE        | <b>Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.</b>               |
|                    | INVALID         | <b>No reportable result. Repeat entire assay run using original stool samples.</b><br>Improper sample preparation, reagent failure, instrument failure or internal control failure.   |
| Negative Control   | POSITIVE        | <b>Incorrect control result. Repeat control testing as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.</b>               |
|                    | NEGATIVE        | Valid negative control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.  |
|                    | INVALID         | <b>No reportable result. Repeat entire assay run using original stool samples.</b><br>Improper sample preparation, reagent failure, instrument failure or internal control failure.   |
| EMPTY WELL         | NONE            | No <i>illumigene</i> Test Device in the <i>illumipro-10</i> Well.<br><b>OR</b><br>The <i>illumigene</i> Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. <b>Repeat the test using original sample.</b> |

## QUALITY CONTROL

- Each device contains an internal control well that controls for amplification inhibition, assay reagents and sample processing effectiveness.
- Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- illumigene C. difficile* External Control Reagents are supplied separately (Catalog 279920). It is recommended that the reactivity of each new lot and each new shipment of *illumigene C. difficile* be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The *illumigene C. difficile* test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate device must be used for each external control reagent.
- DO NOT extract the Positive or Negative Control samples.
- An Extraction Control (a known positive or clinical sample or negative sample spiked with toxigenic *C. difficile*) should be included in each extraction run. The Extraction Control should be treated as a sample.

## EXPECTED VALUES

The frequency of antibiotic-associated diarrhea caused by *C. difficile* is dependent on several factors, including patient population, type of institution and epidemiology. The incidence of CDAD in patients suspected of having antibiotic-associated disease is 15-20% although different facilities may find positive rates above or below this range.<sup>6</sup> The incidence of toxigenic *C. difficile* during the 2010 period of this study was 15.5%.

## LIMITATIONS OF THE PROCEDURE

- Assay does not distinguish between viable and non-viable organisms.
- Test results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures.
- This test is for use with unformed stool samples, preserved or unpreserved, only. Performance characteristic of other clinical specimens types have not been established.
- The performance characteristics for specimens from patients less than two years of age have not been established.

- Two distinct groups have been identified that can harbor *C. difficile* asymptomatically at very high rates. Colonization at rates up to 50% and higher have been reported in infants and rates up to 32% in cystic fibrosis patients.<sup>7,8</sup>
- This assay does not identify antimicrobial susceptibility.
- This test detects but does not differentiate the NAP1 (Ribotype 027) strain from other toxigenic strains of *C. difficile*.
- The detection of bacterial nucleic acid is dependent on proper specimen collection, handling (including transportation and storage) and preparation (dilution and extraction). Failure to follow instructions for collection, handling and preparation may cause incorrect results.

## SPECIFIC PERFORMANCE CHARACTERISTICS

*illumigene C. difficile* was evaluated in 2010 by four independent clinical test sites located in the Midwestern and Southern regions of the United States and the manufacturer. Overall performance information is shown in Table 1.

Table 1. Overall performance data

| Cytotoxic bacterial culture | <i>illumigene C. difficile</i> |              |               |
|-----------------------------|--------------------------------|--------------|---------------|
|                             | Positive                       | Negative     | Total         |
| Positive                    | 99                             | 5**          | 104           |
| Negative                    | 27*                            | 546          | 573           |
| Total                       | 126                            | 551          | 677           |
|                             |                                |              | <b>95% CI</b> |
| <b>Sensitivity</b>          | 99/104                         | <b>95.2%</b> | 89.2 - 97.9%  |
| <b>Specificity</b>          | 546/573                        | <b>95.3%</b> | 93.2 - 96.7%  |
| <b>Correlation</b>          | 645/677                        | <b>95.3%</b> | 93.4 - 96.6%  |

\* 15/27 false positive results were positive by another FDA cleared molecular assay. Of the remaining 12 false positive results, 8 were positive by a FDA cleared assay for the detection of GDH.

\*\* 2/5 false negative results were negative by another FDA cleared molecular assay

A total of 697 qualified patient samples were evaluated. Samples were collected from 274 (39.3%) males and 419 (60.1%) females. In the case of 4 (0.6%) of the patients, sex was not known. The age groups of patients range from 2 years of age to 96 years of age. No differences in test performance were observed based on patient age or sex. Tables 2 and 4 show the assay performance by clinical site and patient age; Table 3 summarizes invalid rate information by site.

Table 2. Performance characteristics by site

| Site   | Positive Samples                                |               |              | Negative Samples                                |               |              |
|--------|---|---------------|--------------|---|---------------|--------------|
|        | <i>illumigene</i> / Cytotoxic bacterial culture | Sensitivity % | 95% CI       | <i>illumigene</i> / Cytotoxic bacterial culture | Specificity % | 95% CI       |
| Total  | 99/104**  | 95.2%         | 89.2 - 97.9% | 546/573*  | 95.3%         | 93.2 - 96.7% |
| Site 1 | 4/5   | 80.0%         | 37.6 - 96.4% | 58/60   | 97.6%         | 88.6 - 99.1% |
| Site 2 | 12/12   | 100%          | 75.7 - 100%  | 62/67   | 92.5%         | 83.7 - 96.8% |
| Site 3 | 20/20   | 100%          | 83.9 - 100%  | 87/92   | 94.6%         | 87.9 - 97.7% |
| Site 4 | 8/8   | 100%          | 67.6 - 100%  | 36/39   | 92.3%         | 79.7 - 97.3% |
| Site 5 | 55/59   | 93.2%         | 83.8 - 97.3% | 303/315   | 96.2%         | 93.5 - 97.8% |

Table 3. Invalid rates by site

| Site   | Clinical Site Evaluation |                  |                     |               |
|--------|--------------------------|------------------|---------------------|---------------|
|        | Total Invalids           | Assay Invalids   | Instrument Invalids | Invalid Rate  |
| Site 1 | 3                        | 3                | 0                   | 3/68 (4.4%)   |
| Site 2 | 1                        | 0                | 1                   | 1/80 (1.3%)   |
| Site 3 | 8                        | 1                | 7                   | 8/120 (6.7%)  |
| Site 4 | 1                        | 1                | 0                   | 1/48 (2.1%)   |
| Site 5 | 7                        | 6                | 1                   | 7/381 (1.8%)  |
| Total  | 20                       | 11/697 (1.6%)*** | 9/697 (1.3%)        | 20/697 (2.9%) |

\*\*\* 1 Specimen remained invalid after repeat testing from the original sample.

Table 4. Results by patient age

| Patient age      | Positive Samples                      |               |              | Negative Samples                      |               |              |
|------------------|---------------------------------------|---------------|--------------|---------------------------------------|---------------|--------------|
|                  | <i>illumigene</i> / Toxigenic culture | Sensitivity % | 95% CI       | <i>illumigene</i> / Toxigenic culture | Specificity % | 95% CI       |
| ≥ 2 - 12 years   | 10/11                                 | 90.9%         | 62.3 - 98.4% | 75/79                                 | 94.9%         | 87.7 - 98.0% |
| > 12 to 21 years | 5/5                                   | 100%          | 56.6 - 100%  | 53/56                                 | 94.6%         | 85.4 - 98.2% |
| > 21 years       | 83/87                                 | 95.4%         | 88.8 - 98.2% | 417/437                               | 95.4%         | 93.0 - 97.0% |
| Age Unknown      | 1/1                                   | 100%          | 20.7 - 100%  | 1/1                                   | 100%          | 20.7 - 100%  |

## ANALYTICAL SENSITIVITY

The analytical sensitivity of this assay for *C. difficile* was based on 20 replicates for each measurand and with a stated probability (eg. 95% where 19/20 replicates are positive) of obtaining positive responses at the following levels of the measurands:

| Strain ID | Toxinotype | Phenotype | LoD/Test    |
|-----------|------------|-----------|-------------|
| VPI 10463 | 0          | A+/B+     | 4 CFU/test  |
| 2007431   | III (NAP1) | A+/B+     | 32 CFU/test |
| CFI       | VIII       | A-/B+     | 64 CFU/test |
| 2006240   | V (NAP7)   | A+/B+     | 32 CFU/test |
| BI8       | III        | A+/B+     | 64 CFU/test |
| 2007858   | IX/XXIII   | A+/B+     | 32 CFU/test |
| 8864      | X          | A-/B+     | 64 CFU/test |

## ASSAY REACTIVITY

The following *C. difficile* stock cultures from different sources were tested and produced positive reactions at 64 CFU/test with *illumigene C. difficile*:

Type 0: Strains 10463, 2004111, 2004205, 2005070, 2005257, 2008029, 2008162, 2008341, 2008351, 2009066, 2009099, B1, G1, J7, K12, Y1; Type III: 2004052, 2004118, 2007431, B17, B18; Type V: 2005325, 2006240, 2008188, 2009018, 2009065, BK6; Type VIII: 43598, 2008016, CF1; Type X 8864; Type XII 2007435; Type IX/XXIII 2007858; Unknown Type 2009132, 2009155, 2009277.

## REPRODUCIBILITY

Blind coded panels of 10 samples were supplied to three independent laboratories for precision studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured at the assay limit of detection (n = 3) and just below the limit of blank (ie, high negative sample, n = 3). The panels also included uncharacterized positive (n = 2) and negative (n = 2) samples. Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of *illumigene C. difficile* were used in this study. The results are given in the table below:

| Sample Type   | Site 1            |      | Site 2            |      | Site 3            |      | Total             |      |
|---------------|-------------------|------|-------------------|------|-------------------|------|-------------------|------|
|               | Percent agreement |      | Percent agreement |      | Percent agreement |      | Percent agreement |      |
| Negative      | 20/20             | 100% | 20/20             | 100% | 19/19****         | 100% | 59/59             | 100% |
| High Negative | 25/30             | 83%  | 29/30             | 97%  | 28/30             | 93%  | 82/90             | 91%  |
| Low Positive  | 30/30             | 100% | 30/30             | 100% | 30/30             | 100% | 90/90             | 100% |
| Positive      | 20/20             | 100% | 20/20             | 100% | 20/20             | 100% | 60/60             | 100% |

\*\*\*\* 1 specimen generated an instrument invalid test result.

## CROSSREACTIVITY STUDIES

Crossreactivity studies were performed with positive and negative stool specimens inoculated with bacterial or fungal organisms to a final concentration of  $1.2 \times 10^8$  or virus at a minimum of  $1 \times 10^{5.06}$  TCID<sub>50</sub>/mL. None of the following organisms in stool reacted with *illumigene C. difficile*:

*Aeromonas hydrophila*, *Bacteroides fragilis*, *Campylobacter coli*, *Campylobacter fetus*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Clostridium sordellii*, *Clostridium perfringens*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* O157:H7, *Escherichia fergusonii*, *Escherichia hermannii*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Listeria monocytogenes*, *Peptostreptococcus anaerobius*, *Plesiomonas shigelloides*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Salmonella* Groups B-E, *Serratia liquefaciens*, *Serratia marcescens*, *Shigella boydii*, *Shigella flexneri*, *Shigella sonnei*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Vibrio parahaemolyticus*, *Yersinia enterocolitica*, Adenovirus Types 40 and 41, Coxsackievirus, Echovirus, Rotavirus.

## TESTS FOR INTERFERING SUBSTANCES

The following substances, at the specified saturated solvent/diluents concentrations, do not interfere with test results in the final concentrations listed: Barium sulfate (5 mg/mL), fecal fat (equivalent to 2.65 mg stearic plus 1.3 mg palmitic acids per mL), hemoglobin (as methemoglobin) (3.2 mg/mL), IgA (5 mg/mL), Imodium AD® (0.00667 mg/mL), Kaopectate® (0.87 mg/mL), Metronidazole (12.5 mg/mL), mucin (3.33 mg/mL) Mylanta® (4.2 mg/mL), Pepto-Bismol® (0.87 mg/mL), Prilosec® (0.5 mg/mL), Tagamet® (0.5 mg/mL), TUMS® (0.5 mg/mL), Vancomycin (12.5 mg/mL), white blood cells (5% V/V), whole blood (5% V/V).

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**REFERENCES**

- Nagamine K, Hase T, Notoni T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. Mol Cell Probes 2002;16:223-29.
- Mori Y, Kitao M, Tomita N, Notoni T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. J Biochem Biophys 2004;59:145-47.
- Rupnik M. Heterogeneity of large clostridial toxins: importance of Clostridium difficile toxinotypes. Microbiol Rev 2008;32:541-55.
- US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
- CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.
- Wilkins TD, Lyerly DM. Clostridium difficile testing: after 20 years, still challenging. J Clin Microbiol 2003;41:531-4.
- Johnson S, Gerding DN. Clostridium-difficile-associated diarrhea. Clin Infect Dis 1998;26:1027-34.
- Lyerly DM, Krivan HC, Wilkins TD. Clostridium difficile: its disease and toxins. Clin Microbiol Rev 1988;1:1-18.

SN11176

REV. 07/10

Manufactured By

Meridian Bioscience, Inc.  
 USA/Corporate Office  
 3471 River Hills Drive  
 Cincinnati, Ohio 45244  
 Telephone: (513) 271-3700  
 Orders/Customer Service:  
 (800) 543-1980  
 Technical Support Center:  
 (800) 343-3858  
 Information Fax: (513) 272-5432  
 Ordering Fax: (513) 271-0124

---

Authorized Representative

Meridian Bioscience Europe  
 Via dell'Industria, 7  
 20020 Villa Cortese (MI)  
 Italy  
 Tel.: +39 0331 433636  
 Fax: +39 0331 433616  
 e-mail: info@mdeur.com

Meridian Bioscience Europe France  
 Le Quadra  
 455, Promenade des Anglais  
 06299 Nice Cedex-3  
 France  
 Tel.: +33 (4) 93 18 72 10  
 Fax: +33 (4) 93 18 72 11  
 e-mail: info@meridianbioscience.fr


Meridian Bioscience Europe s.a. / n.v.  
 Rue de l'Industrie 7  
 1400 Nivelles  
 Belgium  
 Tel.: +32 (0) 67 89 59 59  
 Fax: +32 (0) 67 89 59 58  
 e-mail: info@mdeur.be

Meridian Bioscience Europe b.v.  
 Halderheiweg 6  
 5282 SN Boxtel  
 The Netherlands  
 Tel.: +31 (411) 621166  
 Fax: +31 (411) 624841  
 e-mail: meridian.info@planet.nl

## INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

**Key guide to symbols (Guida ai simboli, Guide des symboles, Guia de símbolos, Erläuterung der graphischen symbole)**

|  |  |  |   |
|--|--|--|---|
|   | Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis  | <b>CONTROL +</b>   | Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle  |
| <b>LOT</b>   | Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung  | <b>CONTROL -</b>   | Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle  |
| <b>IVD</b>   | In vitro diagnostic medical device / Dispositivo médico-diagnóstico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum  | <b>SMP   PREP   DIL   SPE</b>  | Sample Preparation Apparatus containing Sample Diluent / Appareato per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenvorbereitung, in dem sich Probenverdünnungspuffer befindet |
| <b>CE</b>  | This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostika 98/79/EG. | <b>EC REP</b>  | Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft   |
| <b>REF</b>   | Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer  |   | Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren  |
|    | Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten   | <b>RoHS</b>  | Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation de substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe   |
|    | Manufacturer / Fabricante / Fabricant / Fabricante / Hersteller  |   | Caution, consult accompanying documents / Attenzione, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten   |
|    | Contains sufficient for <math>n</math>- tests / Contenuto sufficiente per <math>n</math> saggi / Contenu suffisant pour <math>n</math> tests / Contenido suficiente para <math>n</math> ensayos / Inhalt ausreichend für <math>n</math>-Prüfungen  | <b>BUF   RXN</b>   | Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer   |
|    | Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung   |   | ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifié Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt  |
| <b>SN</b>  | Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer   |   | Recycle – do not dispose of as general waste / Riciclare – non eliminare come rifiuto generico / Recycler – ne pas jeter dans une poubelle / Recycle – no desecha como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen   |
| <b>TEST</b>  | Test Device / Dispositivo test / Dispositif de test / Dispositivo de Proeba / Testgerät  | <b>EX   TUBE</b>   | Extraction tube / Provette per Estrazione / Tubes d'Extraction / Tubo de Extracción / Rohrchen zur Probenextraktion   |
|    | Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum   |   | For IVD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Résultats IVD réservés à l'évaluation des performances / Sólo para evaluación del funcionamiento / Nur zur Leistungsewertung   |
|   | LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Éviter toute exposition au faisceau / Radiación Laser: Evite exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden   |  | HOT SURFACE: Keep hands away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden                                   |
|  | CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement Laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung   | <b>IPX-0</b>   | CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: Vor Feuchtigkeit schützen   |
|  | CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risques de danger / Precaución: Peligroso / WARNUNG: Risikogefahr  |  |   |

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.