PARA-SELLES / KOP-COLOR II FUMOUCZE®

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Nom. : 2000000 - 09/08

PARA-SELLES / KOP-COLOR II FUMOUCZE®

STAINING AND CONCENTRATION OF PARASITIC ELEMENTS IN STOOLS

INTENDED USE:

PARA-SELLES / KOP-COLOR II FUMOUCZE® allows performing:

- the direct examination of parasitic elements in stools with the KOP-COLOR II.
- the concentration of parasitic elements based on Bailenger's method with KOP-COLOR II staining.
- the concentration of parasitic elements based on IODESINE concentration method.

The IODESINE concentration technique replaces the M.I.F. Concentration one. KOP-COLOR II replaces KOP-COLOR. The used staining agents are the same. Base for IODESINE and KOP-COLOR II contain neither formalin nor mercury derivative so as to improve work station safety.

KIT CONTENT:

- | | |
|--------------------------------------|--------------------------|
| - pH 5 aceto-acetate buffer solution | - Conical tubes of 30 mL |
| - Base for IODESINE | - Conical tubes of 10 mL |
| - KOP-COLOR II | - Spatulas |
| - Lugol | - Package insert |

MATERIAL REQUIRED BUT NOT PROVIDED:

- | | |
|--|--------------------------|
| - Micropipette 10-50 µL | - Microscope |
| - Micropipette 100-1000 µL | - Vortex |
| - Haemolysis tubes | - Centrifuge |
| - Pasteur pipettes | - Physiological saline |
| - Slides + coverglasses for microscopy | - Ether or ethyl acetate |
| - Container for contaminated wastes | |

STORAGE CONDITIONS:

Reagents are ready-to-use. Store at +18°...+25°C, **sheltered from sunlight**, until the expiry date indicated on the box. Do not freeze.

SAMPLES STORAGE:

Due to the fragility of some of parasitic stages as protozoa vegetative forms, it is recommended to treat stools as soon as possible after their collection.

WARNINGS AND PRECAUTIONS:

- For in vitro diagnostic use.
- Only for professional use.
- Follow the instructions for use.
- In case of accidental spill of reagents, clean the surface with absorbent paper and water. In case of spill of sample, clean the surface with absorbent paper and bleach.
- Avoid contact of reagent with skin, eyes and mucous membranes. Do not ingest.
- The samples, reagents as well as contaminated materials and products must be eliminated in a container for contaminated wastes, according to the prevailing recommendations and regulations.
- **Base for IODESINE:** R 10: Flammable.

A) DIRECT EXAMINATION AFTER STAINING BY KOP-COLOR II

PRINCIPLE:

KOP-COLOR II is a differential staining process of parasitic elements using a mixture of staining agents one of which is Lugol. Its utilization facilitates the detection of parasitic elements which appear to be yellow, yellow-orange or brownish-yellow on a more or less dark blue background.

TEST PROCEDURE:

- Homogenize stools.
- Take out a volume of stools equivalent to the size of a pea and place it in a haemolysis tube containing 1 mL of thinner (physiological saline, distilled water or pH5 aceto-acetate buffer solution).
- Triturate and shake it to obtain a homogeneous suspension (Vortex shaker).
- Using the micropipette, place 10 µL of KOP-COLOR II on a slide.
- Using the Pasteur pipette, add 1 drop (or, using the micropipette, add 25 µL) of the stools suspension to examine.
- Mix well.
- Place a coverglass over the stools suspension and examine by using a microscope having a white light (blue filter).

INTERPRETATION OF THE RESULTS:

Parasitic elements appear to be yellow, yellow-orange or brownish-yellow on a more or less dark blue background.

B) CONCENTRATION OF PARASITIC ELEMENTS BASED ON BAILENGER'S METHOD WITH KOP-COLOR II STAINING

PRINCIPLE:

Two phase method of concentration using ether as an organic solvent and a pH5 aceto-acetate buffer solution as the aqueous phase. The

examination of sediment is performed after staining by KOP-COLOR II which facilitates the detection of parasitic elements which appear to be yellow, yellow-orange or brownish-yellow on a more or less dark blue background.

TEST PROCEDURE:

- Pour 20 mL of the aceto-acetate buffer in a 30 mL conical tube.
- Homogenize stools.
- Take out the equivalent of a knob of the stools (3 to 4 g or 3 to 4 mL if the stools are fluid) and place it in the aceto-acetate buffer.
- Triturate the stools by using a spatula and shake the mixture vigorously until a homogeneous suspension is obtained (Vortex shaker).
- Permit the mixture to settle for 2 or 3 minutes for the sedimentation of the coarser particles.
- Pour 5 mL of supernatant in a 10 mL conical tube.
- Add 2.5 to 3 mL of ether.
- Plug the tube and shake vigorously in order to obtain an emulsion (shake manually or with a Vortex shaker).
- Take the plug out of the tube and centrifuge the mixture (150-200 g) for 5 minutes to "break" the emulsion.
- In the case of gelling in the higher phase (due to lipophile residues), loosen the emulsion from the tube walls using the Pasteur pipette.
- Eliminate the supernatant by turning the tube upside down.
- Put the sediment in a suspension using 1 or 2 drops of physiological saline (do not allow the sediment to dry up).
- Using the micropipette, place 10 µL of KOP-COLOR II on a slide.
- Using the Pasteur pipette, add 1 drop (or, using the micropipette, add 25 µL) of the suspension to examine.
- Mix well.
- Place a coverglass over the sediment suspension and examine under a microscope using a white light (blue filter).

INTERPRETATION OF THE RESULTS:

Parasitic elements appear to be yellow, yellow-orange or brownish-yellow on a more or less dark blue background.

C) CONCENTRATION OF PARASITIC ELEMENTS BASED ON IODESINE CONCENTRATION METHOD

PRINCIPLE:

Two phase method of concentration using ether as an organic solvent and the IODESINE solution as the aqueous phase. Parasitic elements appear to be pink or more or less darkish brown.

TEST PROCEDURE:

- IODESINE solution preparation:**
 - Place 200 µL of Lugol in the bottom of a 30 mL conical tube, then add 20 mL of Base for IODESINE.
- Homogenize the stools.
- Take out the equivalent of a knob of the stools (3 to 4 g or 3 to 4 mL if the stools are fluid) and place it in the previously prepared Iodesine solution.
- Triturate the stools by using a spatula and shake the mixture vigorously until a homogeneous suspension is obtained (Vortex shaker).
- Permit the mixture to settle for 2 or 3 minutes for the sedimentation of the coarser particles.
- Pour 5 mL of supernatant in a 10 mL conical tube.
- Add 2.5 to 3 mL of ether.
- Plug the tube and shake vigorously to obtain an emulsion (shake manually or with a Vortex shaker).
- Take the plug out of the tube and centrifuge the mixture (150-200 g) for 5 minutes to "break" the emulsion.
- In the case of gelling in the higher phase (due to lipophile residues), loosen the emulsion from the tube walls using the Pasteur pipette.
- Eliminate the supernatant by turning the tube upside down.
- Put the sediment in a suspension using 1 or 2 drops of physiological saline (do not allow the sediment to dry up).
- Using the Pasteur pipette, place 1 drop (or, using the micropipette, add 25 µL) of the suspension to examine on a slide.
- Place a coverglass over the sediment suspension and examine under a microscope using a white light (blue filter). When parasites colouring is too important, it is recommended to increase the white light intensity for a good microscopic observation of these last ones.

INTERPRETATION OF THE RESULTS:

The parasitic elements appear to be pink or more or less darkish brown.

N.B.: POSSIBILITY OF USING ETHYL ACETATE AS ORGANIC SOLVENT INSTEAD OF ETHER

Ether (diethyl ether – C₂H₅ – O – C₂H₅) is the reference solvent for Bailenger and IODESINE concentration methods. It can be replaced by ethyl acetate in equal volume. It is however necessary to note that the solubilisation of certain faecal fragments is more difficult with ethyl acetate than with diethyl ether. **So, with certain stools treated with ethyl acetate, the interface layer may be thicker and the more voluminous sediment may present numerous non-parasitic elements which may hamper microscopic observation.** It is then advisable to take back the sediment with more diluent, for example with 4 drops of physiological water instead of 2, and to make several assemblies for the microscopic examination.

BIBLIOGRAPHIE / REFERENCES:

1. A. O'FEL - Parasitologie mycologie - *Format Utile*, Saint-Maur.
2. J. BAILENGER - Coprologie parasitaire et fonctionnelle - *Imprimerie Drouillard*, Bordeaux.
3. P. BOURÉE - Aide mémoire de parasitologie - *Flammarion*, Paris.
4. A.-M. DELUOL - Atlas de parasitologie - Guide pratique du diagnostic au microscope - (tomes I, II, III), *Edition Varia*, Paris.
5. J.-P. NOZAIS, A. DATRY, M. DANIS, C. BOUDON - Traité de parasitologie médicale - *Pradel*, Paris.
6. M. GENTILLINI, B. DUFLO - Médecine tropicale de voyage - *Flammarion Médecine Sciences*, Paris.
7. Y.-J. GOLVAN - Eléments de parasitologie médicale - *Flammarion*, Paris.
8. H. LEGER, M.-J. NOTTEGHEM - Guide de parasitologie pratique - *SEDES*, Paris.
9. C. JUNOD - Recherche spéciale des oeufs et larves d'Helminthes dans les selles par la méthode des concentrations combinées - *Feuillets de biologie*, 92: 55-62 (1976).
10. D. ENGELS, S. NAHIMANA, B. GRYSELS - Comparison of the direct faecal smear and two thick smear techniques for the diagnosis of intestinal parasitic infections - *Transactions of the Royal Society of Tropical Medicine and Hygiene*, 90: 523-525 (1996).