Mycoplasma DNA Amplification Assay

DNA Amplification Assay for the Detection of *Mycoplasma pneumoniae* in human throat and nasopharyngeal swab specimens

REF 280550

IVD In vitro diagnostic medical device

INTENDED USE

The *illumigene* Mycoplasma DNA amplification assay, performed on the *illumipro-10™*, is a qualitative in vitro diagnostic test for the direct detection of DNA from *Mycoplasma pneumoniae* in human throat and nasopharyngeal swabs obtained from patients suspected of having *Mycoplasma pneumoniae* infection.

The *illumigene* Mycoplasma assay utilizes loop-mediated isothermal DNA amplification (LAMP) technology to detect *Mycoplasma pneumoniae* by targeting a segment of the *Mycoplasma pneumoniae* genome.

Results from the *illumigene* Mycoplasma DNA amplification assay should be used in conjunction with clinical presentation, other laboratory findings, and epidemiological risk factors as an aid in the diagnosis of Mycoplasma infection and should not be used as the sole basis for treatment or other patient management. Positive results do not rule out co-infection with other organisms and negative results in persons with respiratory tract infections may be due to pathogens not detected by this assay. Lower respiratory tract infections due to *M. pneumoniae* may not be detected by this assay. If lower respiratory tract infection due to *M. pneumoniae* is suspected, additional laboratory testing using methods other than the *illumigene* Mycoplasma DNA Amplification Assay may be necessary.

illumigene Mycoplasma is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The *illumigene* Mycoplasma DNA Amplification Assay is based on loop-mediated amplification (LAMP) technology.¹ The assay targets a 208 base pair (bp) sequence of the *Mycoplasma pneumoniae* genome. The target DNA sequence is found in the intracellular protease-like protein gene and is found in all sequences comprising three genomes of *Mycoplasma pneumoniae*.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is the formation of magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian *illumipro-10* Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The *illumigene* Mycoplasma kit includes *illumigene* Assay Control II and *illumigene* Mycoplasma Test Devices. The *illumigene* Assay Control II, used for specimen dilution and preparation, is a buffered solution containing *Staphylococcus aureus* DNA. The *illumigene* Mycoplasma Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Mycoplasma-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA from the Control Reagent and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is diluted with the Control Reagent and combined with *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient specimen allows for parallel processing of target DNA and Control DNA through sample preparation, amplification and detection. The Internal Control monitors DNA extraction, amplification inhibition, assay reagent performance and sample processing effectiveness. The control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The *illumipro-10* monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_f). The *illumipro-10* calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Mycoplasma pneumoniae is a common cause of human upper and lower respiratory infections, including pharyngitis, acute bronchitis and pneumonia.² Mycoplasmas became recognized as human pathogens in the 1960s, with *Mycoplasma pneumoniae* the best known and most studied.³ Characteristics of the organism include a small genome consisting of a single circular chromosome and a triple-layered cell membrane rather than a structured cell wall.³ The size and structure of *Mycoplasma pneumoniae* prevent positive identification by traditional gram stain methods and visualization by light microscopy.

M. pneumoniae has been associated with up to 40%^{2, 3} of community-acquired pneumonias. Infection occurs in both children and adults without geographical, gender or climate-related restrictions.^{3, 6} *M. pneumoniae* is most often associated with atypical pneumonia, presenting with symptoms that include headache, malaise, myalgias, fever, and sore throat accompanied by dry, paroxysmal cough.^{3, 4, 6} While the clinical course of mycoplasmal pneumonia is generally mild and self-limiting, it has been associated with a mortality rate of approximately 1.4%.⁴ An estimated 2 million cases of *M. pneumoniae* infection occur annually with approximately 100,000 pneumonia-related hospitalizations in the United States each year.⁵ It is estimated as many as 18% of children requiring hospitalization with community-acquired pneumonia are caused by *M. pneumoniae*.³

Transmission of *M. pneumoniae* is generally from person to person by aerosol with a reported incubation period of one to four weeks.⁵ As patients with active mycoplasmal infection carry Mycoplasma in the nose, throat, trachea, and sputum, the spread of the disease is facilitated by its accompanying cough.³ This mode of transmission lends itself to outbreaks in close personal contact settings such as schools, military barracks, businesses, summer camps, colleges or institutions. Similarly, the disease is commonly spread among family members within a household.^{3, 4, 5, 6}

Diagnosis of acute *M. pneumoniae* infection is difficult. Basic diagnostic strategies in clinical practice include culture and serology. *M. pneumoniae* culture is often impractical for patient management as the organism may take as many as six weeks to culture.¹ While serology products are widely available and offer practical diagnostic solutions, they may be limited in their ability to detect acute infection. Serological assays are often directed at specific immunoglobulin, IgM and IgG, antibodies in sera. IgM antibodies are generally not detectable with the first seven days of symptom onset and may persist for months after active infection.² Studies indicate that less than 50% of patients with acute infection demonstrate a positive IgG response.²

New diagnostic techniques, including DNA-amplification methods, may enable more rapid diagnosis enabling earlier identification of outbreaks and prevention of secondary cases through implementation of control measures.⁵

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

1. ***illumigene* Assay Control II:** Phosphate-buffered solution containing formalin-treated *E. coli* cells harboring plasmid containing a segment of the *S. aureus* genome with sodium azide (0.09%) as a preservative.
2. ***illumigene* Reaction Buffer II:** Tris-buffered solution containing sodium azide (0.09%) as a preservative.
3. ***illumigene* Mycoplasma Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either Mycoplasma-specific primers (TEST Chamber) or control primers (CONTROL Chamber).
4. **Screw-top Microcentrifuge Tubes:** 1.5 mL RNase/DNase-free microcentrifuge tubes with screw-on caps.

MATERIALS PROVIDED SEPARATELY

1. ***illumigene* Mycoplasma External Control Kit,** Catalog Number: 279940

MATERIALS NOT PROVIDED

1. Qiagen QIAamp® DSP DNA Mini Kit, Catalog Number: 61304
2. Ethanol, 200 Proof
3. Disposable latex gloves, powder free
4. DNase/RNase-free, aerosol resistant pipette tips
5. Specimen collection and transport system
6. Swabs, Nasopharyngeal or Throat: Cotton, Foam, Flocked Nylon, Polyester or Rayon
7. Transport Medium: 0.85% Saline, M4, M4-RT, M5 or UTM-RT

EQUIPMENT NOT PROVIDED

1. Dry-bath with heat block capable of 56 C
2. Dry-bath with 12 mm heat block capable of 95 C
3. Digital thermometer with max/min temperature memory (e.g., Traceable® Lollipop™ waterproof/shockproof thermometer)
4. Vortex mixer
5. Interval timer
6. Microcentrifuge
7. Micropipette(s) capable of dispensing 20, 50, 75, 100, 150, 200 and 500 µL
8. *illumipro-10™*, Meridian Bioscience, Inc. Catalog Number: 610172

PRECAUTIONS

1. All reagents are for in vitro diagnostic use only.
2. Do not interchange Assay Control II Reagent or Test Devices between lots. Reaction Buffer II and Screw-top Microcentrifuge Tubes are interchangeable provided they are within assigned expiration dates when used.

- Follow Biosafety Level 2 and Good Laboratory practices during testing.⁶ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
- Wear disposable gloves while handling specimens and thoroughly wash hands afterwards. Change gloves often.
- Quality Control Programs for Molecular Testing Laboratories should be employed.⁷ To avoid risk of false positives due to organism/amplicon contamination, employ prevention measures (eg, clean work area with suitable cleaner that destroys nucleic acid, use only aerosol resistant pipette tips, handle one specimen at a time) as required by Good Molecular Practices.
- The *illumigene* Mycoplasma Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
- The *illumigene* Mycoplasma Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
- Dispose of used *illumigene* Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

RISK AND SAFETY PHRASES

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-8 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (19-29 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

Refer to the Qiagen QIAamp® DSP DNA Mini Kit package insert for the proper Storage and Handling of Qiagen Mini columns and reagents.

SPECIMEN COLLECTION AND PREPARATION

Sample types: Throat and nasopharyngeal swab specimens.

Sample Collection: Sample collection should be performed in accordance with institutional guidelines for collection of clinical respiratory specimens. Swab sample(s) should be collected with suitable swab types (e.g., Cotton, Foam, Flocked Nylon, Polyester or Rayon).

Place swab(s) in suitable transport medium (e.g., 0.85% Saline, M4, M4-RT, M5 or UTM-RT) and transport to the laboratory. Samples should be held at 2-27 C during transport.

Samples may be held at 19-29 C for up to 48 hours prior to testing. When testing will not be initiated within this time, the sample may be stored at 2-8 C for up to five days.

SPECIMEN PREPARATION

NOTE: Ensure that the *illumipro-10* is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

Specimen Preparation:

- Mix Sample thoroughly.
- Add 50 µL of Assay Control II to an appropriately labeled QIAamp Lysis Tube (LT).
- Add 150 µL of Sample to the Lysis Tube containing Assay Control II.

Specimen Extraction using Qiagen QIAamp® DSP DNA Mini Kit:

NOTE: The steps listed below should be followed for specimen extraction.

- Add 20 µL of QIAamp Proteinase K (PK) to the Lysis Tube containing Assay Control II and Sample.
- Add 200 µL of QIAamp Lysis Buffer (AL) to the Lysis Tube containing Assay Control II, Sample and Proteinase K. Vortex for a minimum of 10 seconds.
- Heat the Lysis Tube at 56 ± 3 C for 10 minutes.
- Add 200 µL Ethanol, 200 Proof, to the Lysis Tube. Vortex for 10 seconds.
- Pipette the entire contents of the Lysis Tube to a QIAamp Mini Spin Column with Wash Tube (WT).
- Centrifuge the Mini Spin Column for 1 minute at ≥ 6000 x g OR 8000 rpm. Discard flow-through and Wash Tube.
- Place the QIAamp Mini Spin Column into a new Wash Tube (WT) and add 500 µL of reconstituted QIAamp Wash Buffer (AW1).
- Centrifuge the Spin Column for 1 minute at ≥ 6000 x g OR 8000 rpm. Discard flow-through and Wash Tube.
- Place the QIAamp Mini Spin Column into a new Wash Tube (WT) and add 500 µL of reconstituted QIAamp Wash Buffer (AW2).
- Centrifuge for 3 minutes at 20,000 x g OR 14,000 rpm.
- Carefully remove the Mini Spin Column from the Wash Tube, taking care to ensure that the column does not come into contact with the flow-through. Discard flow-through and Wash Tube.
- Place the QIAamp Mini Spin Column into an appropriately labeled QIAamp Elution Tube (ET) and pipette 100 µL of QIAamp Elution Buffer (AE) and add to the QIAamp Mini Spin Column without touching the column membrane.
- Hold the QIAamp Mini Spin Column at 19-29 C for 1 minute. Centrifuge for 1 minute at ≥ 6000 x g OR 8000 rpm to elute purified DNA.
- Carefully discard the QIAamp Mini Spin Column.
- Extracted samples may be held at 19-29 C for up to 60 minutes prior to Heat Treatment.

Heat Treatment:

- Heat each Elution Tube containing extracted Specimen/Control DNA in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
- Remove each Elution Tube from the dry-bath/heat block and centrifuge for 10 seconds at 2000 x g OR 500 rpm. Heat-treated samples may be held at 19-29 C for up to 15 minutes prior to testing.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single *illumipro-10* run. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

- Transfer 75 µL of *illumigene* Reaction Buffer II to an appropriately labeled Screw-top Tube (ST Tube).
- Add 75 µL of heat-treated DNA to the Screw-top Tube containing Reaction Buffer II. Vortex for approximately 10 seconds.
- Repeat steps 1 and 2 for all the samples to be tested.
- Remove 1 *illumigene* Mycoplasma Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
- Transfer 50 µL of the heat-treated sample to the TEST chamber (White Bead) of the *illumigene* Test Device. Take care to not introduce extraneous air to the reaction mixture. Using a new pipette tip, transfer 50 µL from the Screw-top Tube containing heat-treated sample to the CONTROL chamber (Yellow Bead) of the *illumigene* Test Device. Do not introduce air bubbles. Close the *illumigene* Test Device and fasten the latch securely.
- Tap device on the bench top or mix to remove air bubbles. Carefully examine the Test Device to ensure that there are no air bubbles left in the tube and no liquid remaining in the top of the device. Amplification and detection should be initiated within 15 minutes.
- Insert the *illumigene* Test Device into the *illumipro-10* and initiate amplification reaction and detection. Please refer to the *illumipro-10* Operator's Manual for detailed instructions. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Mycoplasma pneumoniae</i> target DNA.
	NEGATIVE	No <i>Mycoplasma pneumoniae</i> DNA detected.
	INVALID	No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	NEGATIVE	Incorrect control result. Patient results not reportable. Repeat the control test as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Patient results not reportable. Repeat the control test as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No <i>illumigene</i> Test Device in the <i>illumipro-10</i> Well. OR The <i>illumigene</i> Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

- Each device contains an internal control that controls for amplification inhibition, assay reagents and sample processing effectiveness. Internal control DNA is present in the Assay Control II Reagent and is processed through all steps of the procedure. Primers for amplification of internal control DNA are present in the Control Chamber of the *illumigene* Test Device.

- The heat-treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 ± 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 ± 2 minutes.
- Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
- illumigene* Mycoplasma External Control Reagents are supplied separately (Catalog 279940). It is recommended that the reactivity of each new lot and each new shipment of *illumigene* Mycoplasma be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The *illumigene* Mycoplasma test kit should not be used in patient testing if the external controls do not produce the correct results.
- A separate Test Device must be used for each external control reagent. Refer to the *illumigene* Mycoplasma External Controls package insert for external quality control testing procedures.

EXPECTED VALUES

Overall incidence of *Mycoplasma pneumoniae* in prospectively collected and tested specimens during the 2012 clinical study was 11.7% (12/103).

LIMITATIONS OF THE PROCEDURE

- This product can be used only with the *illumipro-10* instrument.
- The *illumigene* Mycoplasma DNA assay is a qualitative assay and does not provide quantitative values or information about organism load.
- This device has not been evaluated for monitoring treatment of *Mycoplasma pneumoniae* infections.
- Respiratory infections can be caused by *Mycoplasma pneumoniae* as well as other pathogens. Positive results do not preclude coinfection with other respiratory pathogens.
- Test performance has not been established for testing specimens other than human nasopharyngeal and throat swab specimens, for immunocompromised individuals or from patients not suspected of infection with *M. pneumoniae*.
- The detection of *M. pneumoniae* nucleic acid is dependent upon proper specimen collection, handling, transportation, storage and preparation, including extraction. Failure to observe proper procedure in any one of these steps can lead to incorrect results.
- Since there were a limited number of positive specimens in the prospective clinical study, performance characteristics were established with both prospective and retrospective clinical specimens.
- Prevalence of *M. pneumoniae* infection will affect the test's predictive value.
- Organism nucleic acids may persist *in vivo* independent of organism viability. Detection of analyte DNA target may not imply causative agents for clinical symptoms. The *illumigene* Mycoplasma assay does not differentiate between carriers and infected individuals.
- Results should be interpreted together with clinical history, epidemiological data, and other information available to the physician. As with all molecular based diagnostic tests, (A) False negative results may occur from the presence of inhibitors, technical error, sample mix-up or low numbers of organisms in the clinical specimen, (B) False positive results may occur from the presence of cross-contamination by target organisms, their nucleic acids or amplified product, and from non-specific signals.
- The effect of interfering substances was evaluated only for those substances listed in the labeling. Interference by substances other than those described in the "Tests for Interfering Substances" section below could lead to erroneous results.
- Blast analysis does not indicate that the primers react with any other organism except *M. pneumoniae*, however cross-reactivity with respiratory organisms other than those listed in the "Crossreactivity" section below could lead to erroneous results.
- There is potential for a false negative result in the presence of high concentrations of *Moraxella catarrhalis*, *Nocardia asteroides*, or Coronavirus. For these three organisms, a false negative result was observed in only one of seven replicates with samples tested near the Limit of Detection during initial testing that were not confirmed with further testing.
- Phenylephrine HCl found in nasal decongestants produced false negative results at concentrations above 0.595 mg/mL during *M. pneumoniae* strain M129 Limit of Detection replicate testing.

SPECIFIC PERFORMANCE CHARACTERISTICS

The *illumigene* Mycoplasma DNA Amplification Assay was evaluated in 2012 by independent clinical test sites representing geographically distinct regions throughout the United States. A total of 334 qualified throat and nasopharyngeal (NP) swab specimens, collected from patients suspected of *Mycoplasma pneumoniae* infection were evaluated with the test device to establish performance characteristics. Specimens were leftover deidentified specimens submitted to the testing laboratories for routine *M. pneumoniae* testing. Specimens included in performance evaluation were prospective (never frozen) and retrospective (frozen prior to *illumigene* testing). The performance of *illumigene* was compared to a Composite Reference Method that included *M. pneumoniae* bacterial culture with identification and a validated real-time PCR assay followed by bi-directional sequencing for positive specimens. Specimens producing positive *Mycoplasma pneumoniae* results from either the bacterial culture or real-time PCR and bi-directional sequencing were considered positive. Specimens negative for both culture and PCR were considered negative. A total of 12 specimens were excluded from the clinical sample population because culture results were inconclusive and PCR negative, making final disposition of patient status unavailable. A total of 103 (30.8%) prospective and 219 retrospective specimens (65.6%) were tested with one initial invalid result (0.30%).

Tables 1-2 summarize performance characteristics. Statistical analysis of Specimen Type performance data was performed with no significant difference between swab types identified.

Data indicates performance is optimal when specimens are collected and tested prospectively.

Table 1. *illumigene* Mycoplasma Assay Performance: Nasopharyngeal Swabs

Specimen Description	Positive Specimens			Negative Specimens			Invalid Results
	<i>illumigene</i> vs. Comparator	% Sensitivity or PPA	95% CI	<i>illumigene</i> vs. Comparator	% Specificity or NPA	95% CI	
Composite Method Comparator							
Prospective	4/4	% Sensitivity 100%	51.0 – 100.0%	48/48	% Specificity 100%	92.6 – 100.0%	0
Retrospective	34/36	PPA 94.4%	81.9 – 98.5%	86/90	NPA 95.6%	89.1 – 98.3%	0
Validated PCR Assay with Bi-directional Sequencing Comparator							
Prospective	4/4	% Sensitivity 100%	51.0 – 100.0%	48/48	% Specificity 100%	92.6 – 100.0%	0
Retrospective	34/36	PPA 94.4%	81.9 – 98.5%	86/90	NPA 95.6%	89.1 – 98.3%	0

Table 2. *illumigene* Mycoplasma Assay Performance: Throat Swabs

Specimen Description	Positive Specimens			Negative Specimens			Invalid Results
	<i>illumigene</i> vs. Comparator	% Sensitivity or PPA	95% CI	<i>illumigene</i> vs. Comparator	% Specificity or NPA	95% CI	
Composite Method Comparator							
Prospective	8/8	% Sensitivity 100%	67.6 – 100.0%	43/43	% Specificity 100%	91.8 – 100.0%	0
Retrospective	22/26 ^a	PPA 84.6%	66.5 – 93.9%	66/67	NPA 98.5%	92.0 – 99.7%	1
Validated PCR Assay with Bi-directional Sequencing Comparator							
Prospective	8/8	% Sensitivity 100%	67.6 – 100.0%	43/43	% Specificity 100%	91.8 – 100.0%	0
Retrospective	21/21 ^b	PPA 100%	84.5 – 100.0%	70/72	NPA 97.2%	90.4 – 99.2%	1

a. Four specimens originally identified by culture as positive were not confirmed by either the *illumigene* assay or the independent PCR method. Results suggest sample degradation during storage.

b. One specimen originally identified by the *illumigene* assay and culture as positive was negative by the independent PCR method. This specimen is classified as a false positive relative to PCR.

Age information was known for 83.5% (269/322) of the patients included in the performance analysis. Seven (2.6%) patients tested were between 0-28 days of age; 38 (14.1%) patients were between 29 days and up to 2 years of age; 139 (51.7%) patients were between 2 and up to 12 years of age; 61 (22.7%) patients were between 12 and up to 18 years of age; and 9 (3.3%) patients were between 18 and up to 21 years of age. The remaining 15 (5.6%) study patients were 21 years or older. No performance differences were noted based on chronological age.

The study population included 90 (27.9%) female patients and 91 (28.3%) male patients. Gender was unknown for 141 (43.8%) of the study participants. In the specimens for which patient gender was known, no performance differences were noted based on gender.

Clinical performance of the *illumigene* Mycoplasma DNA Amplification Assay was assessed by the testing of deidentified nasopharyngeal and throat swab specimens lacking clinical information; accordingly, the number of patients with *M. pneumoniae* pneumonia included in the clinical studies is unknown and performance for this group cannot be described separately.

ANALYTICAL SENSITIVITY

The analytical sensitivity or Limit of Detection for the *illumigene* Mycoplasma assay was determined for two *M. pneumoniae* strains.

Limit of Detect was determined using a minimum of 20 replicates for each measurand and a stated probability (e.g., 95% where 19/20 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below:

<i>Mycoplasma pneumoniae</i> Strain Description	CFU/Test	CFU/mL
FH (ATCC 15531)	88	2350
M129	7.5	200

ASSAY REACTIVITY

The following *M. pneumoniae* strains were tested and produced positive reactions at or below stated assay limit of detect of 88 CFU/Test (2350 CFU/mL) with *illumigene* Mycoplasma:

PI 1428 (ATCC 29085), MAC (ATCC 15492), M52 (ATCC 15293), Bru (ATCC 15377), M129-B170 (ATCC 29343), Mutant 22 (ATCC 39505) UAB 55612, UAB 56317, UMTB-10G (ATCC 49899). Testing demonstrated that both Type 1 and Type 2 strains are detected with the assay.

REPRODUCIBILITY

Blind-coded panels of 10 samples were supplied to three independent laboratories for reproducibility studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as low positive samples (i.e., near the assay limit of detection, n=3) and high negative samples, (n=3). The panels also included contrived positive (n=3) samples and natural negative samples (n=1). Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of *illumigene* Mycoplasma and five *illumipro-10* instruments were used in this study. Positive and Negative Controls were tested each day of testing. The results are given in the table below:

Sample Type	Site 1		Site 2		Site 3		Total	
	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	
Negative	10/10	100%	10/10	100%	10/10	100%	30/30	100%
High Negative	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Low Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Positive	29/30	96.7%	30/30	100%	29/30	96.7%	88/90	97.8%
Negative Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%
Positive Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%

CROSSREACTIVITY / MICROBIAL INTERFERENCE STUDIES

Crossreactivity studies were performed with positive and negative specimens inoculated with potentially interfering bacterial or fungal organisms to minimum final concentrations of 1.0×10^6 CFU/mL. Viruses were tested at concentrations greater than 1.0×10^5 TCID₅₀/mL or 1.0×10^5 copies/mL. Human DNA was tested at 2.0 ng/test with no crossreactivity observed. Positive specimens contained *M. pneumoniae* concentrations near the limit of detection. None of the following organisms or materials were identified as crossreactive or interfering in the *illumigene* Mycoplasma assay:

Acinetobacter baumannii, *Acinetobacter calcoaceticus*, *Actinomyces odontolyticus*, *Bacillus subtilis*, *Bacteroides fragilis*, *Bordetella parapertussis*, *Bordetella pertussis*, *Burkholderia cepacia*, *Candida albicans*, *Candida glabrata*, *Candida parapsilosis*, *Chlamydia pneumoniae*, *Citrobacter freundii*, *Clostridium difficile*, *Corynebacterium diphtheriae*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Escherichia coli* (ESBL), *Fusobacterium nucleatum*, *Haemophilus ducreyi*, *Haemophilus influenzae*, *Haemophilus parainfluenzae*, *Helicobacter pylori*, *Klebsiella pneumoniae*, *Klebsiella pneumoniae* (KPC), *Legionella pneumophila*, *Listeria monocytogenes*, *Mycoplasma genitalium*, *Mycoplasma hominis*, *Neisseria cinerea*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Peptostreptococcus anaerobius*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella paratyphi* (Group A), *Salmonella typhimurium* (Group B), *Serratia liquefaciens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus agalactiae* (Group B), *Streptococcus anginosus* (Group F), *Streptococcus bovis* (Group D), *Streptococcus canis* (Group G), *Streptococcus equisimilis*, *Streptococcus mitis*, *Streptococcus pneumoniae*, *Streptococcus pyogenes*, *Streptococcus salivarius*, *Ureaplasma urealyticum*, Adenovirus, Coxsackievirus, Cytomegalovirus, Epstein Barr virus, Herpes simplex virus 1, Herpes simplex virus 2, Human metapneumovirus, Influenza A, Influenza B, Parainfluenza virus 1, Parainfluenza virus 2, Parainfluenza virus 3, Respiratory syncytial virus A, Respiratory syncytial virus B, Rhinovirus, Human DNA.

Moraxella catarrhalis, *Nocardia asteroides* and Coronavirus produced unexpected results during initial testing that were not confirmed with further testing. For each organism, false-negative results occurred in one of seven replicates tested near the Limit of Detection.

TESTS FOR INTERFERING SUBSTANCES

Potentially interfering substances were tested with simulated negative and contrived positive (*M. pneumoniae* strains M129 and FH) samples. Potentially interfering substances were diluted in sterile saline and added to M4 medium with rayon swabs (Negative Sample) and to M4 medium with polyester swabs (Contrived Positive Sample and tested).

The following biological substances, at the saturated solvent/diluent concentrations indicated, do not interfere with *illumigene* Mycoplasma test results:
Mucus (5.0mg/mL), White Blood Cells (0.5% v/v) and Whole Blood (5% v/v).

The following chemical substances, at the saturated solvent/diluents concentrations indicated, do not interfere with *illumigene* Mycoplasma test results:

Acetaminophen (18.1 mg/mL), Albuterol Sulfate (20 mg/mL), Aspirin (9.1 mg/mL), Azithromycin dehydrate (2.0 mg/mL), Cepacol® Mouthwash [Ethanol, denatured (1.4% v/v), Cetylpyridinium chloride (0.005% v/v)], Contac® Cold + Flu Tablets [Acetaminophen (14.8 mg/mL), Chlorpheniramine maleate (0.06 mg/mL), Phenylephrine HCl (0.15 mg/mL)], Diphenhydramine HCl (2.6 mg/mL), Erythromycin (20.0 mg/mL), HALLS® Cough Drops [Menthol (0.06 mg/mL)], Ibuprofen (12.7 mg/mL), Phenylephrine HCl (0.595 mg/mL), Prednisone (20.0 mg/mL), Robitussin® Cough+Chest Congestion Syrup [Dextromethorphan HBr (0.20 mg/mL), Guaifenesin (2.0 mg/mL)], Saline Nasal Spray [Sodium chloride (0.65 mg/mL)].

Phenylephrine HCl found in nasal decongestants produced false negative results at concentrations above 0.595 mg/mL during *M. pneumoniae* strain M129 Limit of Detection replicate testing.

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REFERENCES

1. Nagamine, K., Hase, T., Notoni, T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002;16:223-29.
2. Waites, K.B., et al. New insights into the pathogenesis and detection of *Mycoplasma pneumoniae* infections. *Future Microbiol.* 2008 December; 3(6):635-648.
3. Waites, K.B., Talkington, D.F., *Mycoplasma pneumoniae* and Its Role as a Human Pathogen. *Clin Microbiol Rev.* 2004 October; 17(4):697-728.
4. Thibodeau, K.P., Viera, A.J. Atypical Pathogens and Challenges in Community-Acquired Pneumonia. *American Family Physician.* 2004 April; Vol 66, Number 7.
5. US Department of Health and Human Services, Centers for Disease Control and Prevention. *Mycoplasma pneumoniae*. 2005; http://www.cdc.gov/ncidod/dbmd/diseaseinfo/mycoplasmapneum_t.htm
6. American Academy of Pediatrics. *Red Book: 2012 Mycoplasma pneumoniae and Other Mycoplasma Species Infections*. Pickering LK, ed. 29th ed. Elk Grove Village, IL: American Academy of Pediatrics; 2012.
7. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.

SN11183

REV. 05/13

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: (513) 271-3700
Orders/Customer Service:
(800) 543-1980
Technical Support Center:
(800) 343-3858
Information Fax:(513) 272-5432
Ordering Fax:(513) 271-0124

Manufactured By

Meridian Bioscience Europe
Via dell'Industria, 7
20020 Villa Cortese (MI)
Italy
Tel: +39 0331 433636
Fax: +39 0331 433616
e-mail: info@mdeur.com

Meridian Bioscience Europe s.a. / n.v.
Rue de l'Industrie 7
1400 Nivelles
Belgium
Tel.: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
e-mail: info@mdeur.be

Meridian Bioscience Europe b.v.
Halderheiweg 6
5282 SN Boxtel
The Netherlands
Tel.: +31 (411) 621166
Fax: +31 (411) 624841
e-mail: meridian.info@planet.nl

Meridian Bioscience Europe France
34, rue de Ponthieu
75008 Paris
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
e-mail: info@meridianbioscience.fr

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product:

Key guide to symbols (Guida ai simboli, Guide des symbols, Guia de símbolos, Erläuterung der graphischen symbole)

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For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.