Group A *Streptococcus* (Group A Strep) DNA Amplification Assay

DNA Amplification Assay for the Detection of *Streptococcus pyogenes* in throat swab samples

REF 280150

IVD In vitro diagnostic medical device

INTENDED USE

The *illumigene* Group A *Streptococcus* (Group A Strep) assay, performed on the *illumipro-10*[™], is a qualitative in vitro diagnostic test for the detection of *Streptococcus pyogenes* (Group A β -hemolytic *Streptococcus*) in throat swab specimens.

The *illumigene* Group A Strep assay utilizes loop-mediated isothermal DNA amplification (LAMP)^{1,2} technology to detect *Streptococcus pyogenes* by targeting a segment of the *Streptococcus pyogenes* genome. Results from the *illumigene* Group A Strep assay can be used as an aid in the diagnosis of Group A Streptococcal pharyngitis. The assay is not intended to monitor treatment for Group A *Streptococcus* infections.

illumigene Group A Strep is intended for use in hospital, reference or state laboratory settings. The device is not intended for point-of-care use.

SUMMARY AND EXPLANATION OF THE TEST

The *illumigene* Group A *Streptococcus* molecular assay is based on loop-mediated amplification (LAMP) technology. The assay targets a highly conserved 206 base pair (bp) sequence of the *Streptococcus pyogenes* genome residing in a region of the pyrogenic exotoxin B (speB) gene.

Loop-mediated amplification uses specially designed primers to provide for specific and continuous isothermal DNA amplification. A by-product of amplification is magnesium pyrophosphate, which forms a white precipitate leading to a turbid reaction solution. Reaction solution absorbance characteristics are monitored by the Meridian *illumipro-10* Incubator/Reader. Changes in reaction solution absorbance characteristics created by precipitation of magnesium pyrophosphate indicate the presence of target DNA. The absence of target DNA results in no significant change in sample absorbance.

The *illumigene* Group A Strep kit includes *illumigene* Sample Preparation Apparatuses and *illumigene* Group A *Streptococcus* Test Devices. The *illumigene* Sample Preparation Apparatus II, used for specimen dilution and preparation, is filled with a buffered solution containing formalin treated *E. coli* harboring *Staphylococcus aureus* DNA. The *illumigene* Group A Strep Test Device contains one lyophilized amplification reagent bead in each of two chambers: a TEST chamber with Group A *Streptococcus*-specific primers and a CONTROL chamber with *S. aureus*-specific primers. The *S. aureus* DNA in the Sample Preparation Apparatus and the *S. aureus*-specific primers in the CONTROL chamber function as the Internal Control for the assay. During specimen preparation, each patient specimen is added to the Sample Preparation Apparatus II and combined with the *S. aureus* DNA prior to amplification. Addition of *S. aureus* DNA to the patient sample allows for parallel processing of target DNA and Control DNA through amplification and detection. The Internal Control monitors amplification inhibition, assay reagent performance and sample processing effectiveness. The Control *S. aureus* target must be amplified and detected in the final reaction or the test is considered invalid and patient results are not reported.

The *illumipro-10* monitors changes in absorbance characteristics by measuring transmission of light through the Test and Control reaction solutions. Light transmission is checked at the assay Run Start (Signal_{initial}, S_i) and at the assay Run End (Signal_{final}, S_f). The *illumipro-10* calculates the change in light transmission between Run End and Run Start (S_f:S_i) and compares the ratio to a fixed cut-off value.

Fixed cut-off values for the TEST chamber are used to report sample results. TEST chamber S_f:S_i ratios less than 82% are reported as 'POSITIVE'; TEST chamber S_f:S_i ratios greater than or equal to 82% are reported as 'NEGATIVE'. Numerical values are not reported.

Fixed cut-off values for the CONTROL chamber are used to determine validity. CONTROL chamber S_f:S_i ratios less than 90% are considered valid and allow for reporting of TEST chamber results (POSITIVE, NEGATIVE). CONTROL chamber S_f:S_i ratios greater than or equal to 90% are considered invalid and prevent reporting of TEST chamber results. Invalid CONTROL chamber reactions are reported as 'INVALID'. Numerical values are not reported.

More stringent cut-off criteria are applied to the CONTROL chamber reaction to ensure amplification is not inhibited, reagents are performing as intended and that sample processing was performed appropriately.

BIOLOGICAL PRINCIPLES

Group A *Streptococcus*, or *Streptococcus pyogenes*, is a bacterium commonly found in the human throat or on skin. *S. pyogenes* causes a wide variety of diseases in humans, the most common being acute pharyngitis or strep throat. Pharyngitis is diagnosed in approximately 11 million patients in the United States each year, with Group A *Streptococcus* accounting for 15% - 30% of cases in children and 5% - 20% of cases in adults.³ Streptococcal pharyngitis is often accompanied by sore throat, fever, tonsillar exudates and enlarged cervical lymph nodes.

Group A *Streptococcus* infections may result in mild illness (e.g. pharyngitis, impetigo) or may lead to invasive, life-threatening illness (e.g. cellulitis, bacteremia, necrotizing fasciitis, Streptococcal toxic shock syndrome). Approximately 25% of patients with necrotizing fasciitis will die from the infection; more than 35% of patients with streptococcal toxic shock syndrome will die from the infection.⁴ *Streptococcus pyogenes* may also present in healthy, asymptomatic patients. Throat culture surveys of healthy children taken during school outbreaks have shown a carrier prevalence of up to 20%.³

Signs and symptoms of Group A *Streptococcus* pharyngitis often overlap with other infections, making laboratory diagnosis an important tool in identification and treatment of infection. Specimens for laboratory testing should be obtained by vigorous swabbing of both tonsils and the posterior pharynx.⁵ Traditional laboratory diagnosis is performed by rapid antigen testing or culture on sheep blood agar followed by differentiation with latex agglutination.⁵ Specificities of rapid antigen tests is variable and confirmation of negative results with follow-up throat swab culture is recommended.⁵ Throat swab cultures may be 90 to 95% sensitive when correct sampling and plating techniques are used.³

Diagnostic tests such chemiluminescent DNA probes may be as sensitive as standard throat cultures on sheep blood agar.³

REAGENTS/MATERIALS PROVIDED

The maximum number of tests obtained from this test kit is listed on the outer box.

- illumigene* Sample Preparation Apparatus II/Negative Control III (SMP PREP):** Tris-buffered solution containing formalin-treated *E. coli* harboring *S. aureus* DNA and sodium azide (0.09%) as a preservative.
- illumigene* Group A Strep Test Device:** Two-chambered device containing lyophilized amplification reagents (DNA polymerase, deoxynucleotide triphosphates) and either Group A *Streptococcus*-specific primers (TEST Chamber) or *Staphylococcus aureus*-specific primers (CONTROL Chamber).
- illumigene* Heat Treatment Tubes**

MATERIALS PROVIDED SEPARATELY

illumigene Group A *Streptococcus* (Group A Strep) External Control Kit, Catalog Number: 279910

MATERIALS NOT PROVIDED

- Disposable latex gloves, powder free
- DNase/RNase-free, aerosol resistant pipette tips
- Specimen collection and transport system
- Swabs, breakable plastic shaft: Rayon, Polyester or Flocked Nylon
- Non-nutritive Transport Medium: Liquid Amies, without charcoal; Liquid Stuart

EQUIPMENT NOT PROVIDED

- Dry-bath with 12 mm heat block capable of 95 C
- Digital thermometer with Max/Min Temperature Memory (eg, Traceable® Lollipop[™] Waterproof/Shockproof Thermometer)
- Vortex mixer
- Interval timer
- Micropipette capable of dispensing 50 μ L
- illumipro-10*[™], Meridian Bioscience, Inc. Catalog Number: 610172

PRECAUTIONS

- All reagents are for in vitro diagnostic use only.
- Do not interchange Sample Preparation Apparatuses or Test Devices between lots. Heat Treatment Tubes are interchangeable provided they are within assigned expiration date when used.
- Follow Biosafety Level 2 and Good Laboratory practices during testing.⁶ Treat all specimens and used Test Devices as capable of transmitting infectious agents. Do not eat, drink or smoke in areas where specimens or kit reagents are handled.
- Wear disposable gloves while handling specimens and thoroughly wash hands afterwards.
- Quality Control Programs for Molecular Testing Laboratories, including proper use and care of equipment, should be employed.⁷
- The *illumigene* Group A Strep Test Device contains lyophilized reagents. The protective pouch should not be opened until ready to perform the assay.
- The *illumigene* Group A Strep Test Device includes a latch feature that is designed to prevent contamination of the test area with amplification product. Do NOT use Test Devices with broken latches.
- Dispose of used *illumigene* Test Devices immediately after processing, leaving the device latch securely in place. Do NOT open the Test Device after processing. Opening the device after amplification may result in contamination of the test area with amplification product.

RISK AND SAFETY PHRASES

There are no known hazards associated with this product.

SHELF LIFE AND STORAGE

The expiration date is indicated on the kit label. Store the kit at 2-27 C.

REAGENT PREPARATION

Ensure kit reagents are at room temperature (21-27 C) before use. Incorrect results may be obtained if reagents are not brought to room temperature prior to use.

SPECIMEN COLLECTION AND PREPARATION

Sample type: Throat swabs

Sample Collection: Throat sample collection should be performed in accordance with institutional guidelines for collection of clinical specimens for culture of Group A *Streptococcus*. Samples should be collected by vigorously swabbing the tonsils and the posterior pharynx.

Place swab(s) in a non-nutritive transport medium (e.g. Liquid Amies, without charcoal or Liquid Stuart) and transport to the laboratory. Samples should be held at 2-27 C during transport.

Samples may be held at 21-27 C for up to 48 hours prior to testing. When testing will not be initiated within this time, the sample may be stored at 2-8 C for up to seven days.

SPECIMEN PREPARATION

NOTE: Ensure that the *illumipro-10* is powered on and required performance verifications have been completed prior to initiation of SPECIMEN PREPARATION. Refer to the *illumipro-10* Operator's Manual for further information regarding instrument set-up and operation.

1. Add the throat swab to the *illumigene* SMP PREP. Break off the handle of the swab. Replace and secure SMP PREP cap. Vortex the SMP PREP for a minimum of 10 seconds. Specimens in the *illumigene* SMP PREP may be held at 2-29 C for up to 2 hours prior to heat treatment.
2. Remove the tip cap from the SMP PREP and squeeze five to 10 drops of sample into a clean *illumigene* Heat Treatment Tube.
3. Repeat Sample Preparation Steps for all samples to be processed.
4. Heat each Sample/Control mixture in a dry-bath/heat block at 95 ± 5 C for 10 ± 2 minutes. Monitor heat-treatment step with digital thermometer and interval timer.
5. Remove each Heat Treatment Tube from the dry-bath/heat. Heat-treated samples may be held at 21-29 C for up to one hour prior to testing.

TEST PROCEDURE

NOTE: A maximum of 10 samples can be processed in a single *illumipro-10* run.

1. Vortex heat-treated samples for approximately 10 seconds.
2. Remove 1 *illumigene* Group A Strep Test Device from its protective pouch per sample. Carefully open the device, holding the chambers such that the lyophilized reagent will not fall out upon opening. Place device on a flat surface or in a rack that can accommodate the device.
3. Transfer 50 µL of the heat-treated sample to the TEST chamber (White Bead) of the *illumigene* Test Device. Take care to not introduce extraneous air to the reaction mixture. Using a new pipette tip, transfer 50 µL of the heat-treated sample to the CONTROL chamber (Yellow Bead) of the *illumigene* Test Device. Take care to not introduce extraneous air to the reaction mixture. Close the *illumigene* Test Device and fasten the latch securely.
4. Tap device on the bench top or mix to remove air bubbles. Carefully examine the test device for dissolution of the Control/Test Bead, for air bubbles left in the tube and liquid in the top of the device. If undissolved beads, air bubbles or liquid is noted, tap the device on the bench top and repeat the visual inspection.
5. Insert the *illumigene* Test Device into the *illumipro-10* and initiate amplification reaction and detection. Results will be displayed at the conclusion of the run.

INTERPRETATION OF RESULTS

Sample ID	Reported Result	Interpretation
Patient Specimen	POSITIVE	Sample contains <i>Streptococcus pyogenes</i> target DNA.
	NEGATIVE	No <i>Streptococcus pyogenes</i> DNA detected.
	INVALID	No reportable result. Repeat the test using the original sample. Inhibitory patient specimen, improper sample preparation, reagent failure, instrument failure or internal control failure.
Positive Control	POSITIVE	Valid positive control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	NEGATIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
Negative Control	POSITIVE	Incorrect control result. Repeat the control tests as the first step in determining the root cause of the failure. If control failures are repeated please contact Meridian's Technical Services at 1-800-343-3858 (US) or your local distributor.
	NEGATIVE	Valid negative control result. Reagents active at time of use, <i>illumipro-10</i> performing correctly.
	INVALID	No reportable result. Repeat entire assay run using original samples. Improper sample preparation, reagent failure, instrument failure or internal control failure.
EMPTY WELL	NONE	No <i>illumigene</i> Test Device in the <i>illumipro-10</i> Well. OR The <i>illumigene</i> Test Device present is compromised due to sample preparation failure, dirty device or improperly seated device. Repeat the test using original sample.

QUALITY CONTROL

This test should be performed per applicable local, state, or federal regulations or accrediting agencies.

1. Each device contains an internal control chamber that controls for amplification inhibition, assay reagents and sample processing effectiveness.
2. The heat treatment step is monitored with an external thermometer and interval timer. Use the max/min temperature memory of the thermometer to ensure that a temperature of 95 ± 5 C is maintained. Use the interval timer to ensure that heat-treatment duration is 10 ± 2 minutes.
3. Good laboratory practice recommends the use of control materials. Users should follow the appropriate federal, state and local guidelines concerning the running of external quality controls.
4. *illumigene* Group A *Streptococcus* External Control Reagents are supplied separately (Catalog 279910). It is recommended that the reactivity of each new lot and each new shipment of *illumigene* Group A Strep be verified on receipt and before use. External control tests should be performed thereafter in accordance with appropriate federal, state and local guidelines. The *illumigene* Group A Strep test kit should not be used in patient testing if the external controls do not produce the correct results.
5. A separate device must be used for each external control reagent.

EXPECTED VALUES

Overall incidence of *Streptococcus pyogenes* in patients tested during the 2012 clinical study was 14.6% (116/796).

LIMITATIONS OF THE PROCEDURE

1. The *illumigene* Group A Strep assay does not distinguish between viable and non-viable organisms.
2. Respiratory infections can be caused by Streptococci of serogroups other than A as well as other pathogens. This device does not differentiate between carriers and infected individuals.
3. The use of antibiotics or over-the-counter medications may suppress the *Streptococcus* growth in culture despite the presence of organisms detectable by nucleic acid tests.
4. As with all diagnostic tests, results should be interpreted together with other information available to the physician.
5. Zincum Aceticum 2X, Zincum Gluonicum 1X, as found in Zicam® Homeopathic Cold Remedies, interfere with the *illumigene* Group A Strep Assay.

SPECIFIC PERFORMANCE CHARACTERISTICS

The *illumigene* Group A *Streptococcus* DNA Amplification Assay was evaluated in 2012 by independent clinical test sites representing geographically distinct regions throughout the United States. Compiled clinical study information is described below.

Table 1. *illumigene* Group A Strep Performance Data Summary

Group A <i>Streptococcus</i> Composite Culture	<i>illumigene</i> Group A Strep		
	Positive	Negative	Total
Positive	100	2	102
Negative	16	680*	696
Total	116	682	798
95% CI			
Sensitivity	100/102	98.0%	93.1 – 99.5%
Specificity	680/696	97.7%	96.3 – 98.6%
Invalid Rate	2/800	0.3%	0.1 – 0.9%

* Two specimens produced initial invalid results that were negative upon re-test. A total of 798 qualified specimens were evaluated with the test device and composite bacterial culture methods.

Composite Culture Method testing included the Clinical Site Culture Method and a Reference Culture Method performed by Meridian Bioscience. Site culture testing was performed by direct plating of the throat swab specimen while reference culture testing was performed by plating the swab transport media. Specimens producing positive Group A *Streptococcus* results from either the Site Culture Method or the Reference Culture Method were considered positive.

Table 2 and Table 3 summarize performance information by Clinical Site. Site data was analyzed based on Composite Culture Methods and the Site Culture Method. Seven of the 102 Composite Culture-positive specimens were positive by Site Culture methods only; 26 of the 102 Composite Culture-positive specimens were positive only by the Reference Culture method; the two *illumigene* false-negative results were positive only by the Reference Culture performed by Meridian Bioscience. Statistical analysis of Site performance data was performed with no significant difference among Sites identified.

Table 2: *illumigene* Group A Strep Assay Performance by Site; Composite Culture Method

	Positive Specimens			Negative Specimens		
	<i>illumigene</i> /Composite Culture	% Sensitivity	95% CI	<i>illumigene</i> /Composite Culture	% Specificity	95% CI
Total	100/102	98.0%	93.1 - 99.5%	680/696	97.7%	96.3 - 98.6%
Site 1	47/47	100.0%	92.4 - 100.0%	287/291	98.6%	96.5 - 99.5%
Site 2	28/30	93.3%	78.7 - 98.2%	204/212	96.2%	92.7 - 98.1%
Site 3	25/25	100.0%	86.7 - 100.0%	189/193	97.9%	94.8 - 99.2%

Table 3: illumigene Group A Strep Assay Performance by Site; Site Culture Method

	Positive Specimens			Negative Specimens		
	illumigene /Site Culture	% Sensitivity	95% CI	illumigene /Site Culture	% Specificity	95% CI
Total	74/74	100.0%	95.1 - 100.0%	682/724	94.2%	92.3 - 95.7%
Site 1	40/40	100.0%	91.2 - 100.0%	287/298	96.3%	93.5 - 97.9%
Site 2	18/18	100.0%	82.4 - 100.0%	206/224	92.0%	87.7 - 94.9%
Site 3	16/16	100.0%	80.6 - 100.0%	189/202	93.6%	89.3 - 96.2%

Seventy two (9.0%) patients tested were two years of age or younger; 385 (48.2%) patients were between two and 12 years of age; while 259 (32.5%) patients were greater than 12 and less than 21 years of age. The remaining 82 (10.3%) study patients were greater than 21 years old. Age information was known for all patients included in the performance analysis. No performance differences were noted based on chronological age.

The study population included 410 (51.4%) female patients and 386 (48.4%) male patients. Gender was unknown for two (0.2%) of the study participants. No performance differences were noted based on gender.

ANALYTICAL SENSITIVITY

The analytical sensitivity or Limit of Detect for the illumigene Group A Strep assay was determined for common *S. pyogenes* strains.

Limit of Detect was determined using a minimum of 20 replicates for each measurand and a stated probability (eg. 95% where 19/20 replicates are positive) of obtaining positive responses. Analytical sensitivity testing is summarized below:

Streptococcus pyogenes Strain Description	CFU/Test
ATCC 12344	400
ATCC 19615	430

ASSAY REACTIVITY

The following *S. pyogenes* strains were tested and produced positive reactions at or below stated assay limit of detect of 400 CFU/Test with illumigene Group A Strep: ATCC12384, NCIMB 13285, CCUG 33061, CCUG 33409, CCUG 39158, ATCC 49399, CCUG 53553.

REPRODUCIBILITY

Blind-coded panels of 10 samples were supplied to three independent laboratories for reproducibility studies. Samples were randomly sorted within each panel to mask sample identities. The panels included contrived samples manufactured as low positive samples (ie. near the assay limit of detection, n=3) and high negative samples, (n=3). The panels also included contrived positive (n=3) samples and natural negative samples (n=1). Testing was performed by different operators at each site on the same day (intra-assay variability) for five days (inter-assay variability). Three lots of illumigene Group A Strep and five illumipro-10 instruments were used in this study. Positive and Negative Controls were tested each day of testing. The results are given in the table below:

Sample Type	Clinical Site 1		Clinical Site 2		Clinical Site 3		Total	
	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	Percent agreement	
Negative	10/10	100%	10/10	100%	10/10	100%	30/30	100%
High Negative	29/30	96.7%	30/30	100%	28/30	93.3%	87/90	96.7%
Low Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Positive	30/30	100%	30/30	100%	30/30	100%	90/90	100%
Negative Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%
Positive Control	10/10	100%	10/10	100%	10/10	100%	30/30	100%

CROSSREACTIVITY

Crossreactivity studies were performed with positive and negative specimens inoculated with bacterial or fungal organisms to a final concentration of 1.2 x 10⁸ CFU/mL. Human DNA was tested at 0.02 mg/mL with no cross-reactivity observed. None of the following organisms reacted with illumigene Group A Strep:

Acinetobacter baumannii, *Acinetobacter lwoffii*, *Aeromonas hydrophila*, *Arcanobacterium haemolyticum*, *Bordetella bronchiseptica*, *Bordetella holmesii*, *Bordetella parapertussis*, *Bordetella pertussis*, *Burkholderia cepacia*, *Campylobacter jejuni*, *Candida albicans*, *Citrobacter freundii*, *Corynebacterium diphtheria*, *Corynebacterium pseudodiphtheriticum*, *Enterococcus faecalis*, *Enterococcus faecium*, *Escherichia coli*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Lactococcus lactis*, *Legionella jordanii*, *Legionella micdadei*, *Legionella pneumophila*, *Listeria monocytogenes*, *Moraxella catarrhalis*, *Morganella morganii*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Pseudomonas fluorescens*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Stenotrophomonas maltophilia*, *Streptococcus agalactiae*, *Streptococcus anginosus*, *Streptococcus bovis*, *Streptococcus canis*, *Streptococcus dysgalactiae* (subspecies *equisimilis*), *Streptococcus equinus*, *Streptococcus intermedius*, *Streptococcus mitis*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus salivarius*, *Streptococcus suis*, *Streptococcus uberis*, *Streptococcus sp. Viridans* type.

Two organisms produced unexpected results during initial testing that were not confirmed through further testing. *Arcanobacterium haemolyticum* produced false-negative results for one of 13 replicates tested. *Haemophilus influenzae* produced false-positive results for one of 13 replicates tested.

Mycoplasma pneumoniae was tested at 1.5 x 10⁶ CFU/mL with no reaction with the illumigene Group A Strep assay.

TESTS FOR INTERFERING SUBSTANCES

Interfering substance testing was performed by adding potentially interfering substances diluted in sterile saline or contrived positive samples (ATCC 12344, ATCC 19615) to the illumigene SMP PREP. Rayon swabs inoculated with negative throat flora and Liquid Amies Transport Media were added to the illumigene SMP PREP containing potentially interfering substances and tested.

The following biological substances, at the saturated solvent/diluent concentrations indicated, do not interfere with test results: Mucus (5.0mg/mL), Human saliva (10% v/v), and Whole Blood (2.5% v/v). Whole Blood at concentrations greater than 2.5% v/v may interfere with the illumigene Group A Strep assay

The following chemical substances, at the saturated solvent/diluents concentrations indicated, do not interfere with test results:

Acetaminophen (19.5 mg/mL), Aspirin (12.3 mg/mL), Cepacol® Mouthwash, [Cetylpyridinium Chloride (0.005% v/v)], Cepacol® Sore Throat Lozenges [Benzocaine (0.09 mg/mL), Menthol (0.02 mg/mL)], Chloraseptic® Oral Anesthetic/Analgesic [Phenol (0.07% v/v)], Contac® Cold + Flu Tablets [Acetaminophen (16.2 mg/mL), Chlorpheniramine maleate (0.06 mg/mL), Phenylephrine HCl (0.16 mg/mL)], Crest® Complete Toothpaste [Sodium fluoride (0.1 mg/mL)], Diphenhydramine HCl (2.7 mg/mL), HALLS® Cough Drops [Menthol (0.08 mg/mL)], Ibuprofen (15.6 mg/mL), Listerine® Antiseptic Mouthwash [Eucalyptol (0.0092% v/v), Menthol (0.0042% v/v), Methyl salicylate (0.0060% v/v), Thymol (0.0064% v/v)], Robitussin® Cough/Chest Congestion Cough Syrup [Dextromethorphan HBr (0.2 mg/mL), Guaifenesin (2.0 mg/mL)].

Zicam® Oral Mist [Zincum Aceticum 2X, Zincum Gluonicum 1X (0.625% v/v)] produced invalid results in all replicates tested.

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REFERENCES

- Nagamine, K., Hase, T., Notoni, T. Accelerated reaction by loop-mediated isothermal amplifications using loop primers. *Mol Cell Probes* 2002;16:223-29.
- Mori, Y., Kitao, M., Tomita, N., Notoni, T. Real-time turbidimetry of LAMP reaction for quantitating template DNA. *J Biochem Biophys* 2004;59:145-47.
- Choby, B.A. Diagnosis and treatment of Streptococcal Pharyngitis. *Am Fam Phys*: 2009; 79(5):383-90.
- US Department of Health and Human Services, Centers for Disease Control and Prevention. Group A Streptococcal Disease. www.cdc.gov; 2008.
- American Academy of Pediatrics. *Red Book: 2012 Report of the Committee on Infectious Diseases*. Pickering LK, ed. 29th ed. Elk Grove Village, IL: American Academy of Pediatrics; 2012. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
- CLSI: MM3-A2 Molecular diagnostic methods for infectious disease; approved guideline, 2nd ed. Wayne PA: Clinical Laboratory Standards Institute. 2006.

SN11187

REV. 09/12

Meridian Bioscience, Inc.
USA/Corporate Office
3471 River Hills Drive
Cincinnati, Ohio 45244
Telephone: (513) 271-3700
Orders/Customer Service:
(800) 543-1980
Technical Support Center:
(800) 343-3858
Information Fax:(513) 272-5432
Ordering Fax:(513) 271-0124

Meridian Bioscience Europe
Via dell'Industria, 7
20020 Villa Cortese (MI)
Italy
Tel: +39 0331 433636
Fax: +39 0331 433616
e-mail: info@mdeur.com

Meridian Bioscience Europe s.a. / n.v.
Rue de l'Industrie 7
1400 Nivelles
Belgium
Tel.: +32 (0) 67 89 59 59
Fax: +32 (0) 67 89 59 58
e-mail: info@mdeur.be


Meridian Bioscience Europe b.v.
Halderheiweg 6
5282 SN Boxtel
The Netherlands
Tel: +31 (411) 621166
Fax: +31 (411) 624841
e-mail: meridian.info@planet.nl

Meridian Bioscience Europe France
34, rue de Ponthieu
75008 Paris
Tel: +33 (0) 1 42 56 04 40
Fax: +33 (0) 9 70 06 62 10
e-mail: info@meridianbioscience.fr

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product.

Key guide to symbols (Guida ai simboli, Guide des symbols, Guia de simbolos, Erläuterung der graphischen symbole)

	Use By / Utilizzare entro / Utiliser jusque / Fecha de caducidad / Verwendbar bis	CONTROL +	Positive control / Controllo positivo / Contrôle positif / Control positivo / Positive Kontrolle
LOT	Batch Code / Codice del lotto / Code du lot / Código de lote / Chargenbezeichnung	CONTROL -	Negative control / Controllo negativo / Contrôle négatif / Control negativo / Negative Kontrolle
IVD	In vitro diagnostic medical device / Dispositivo medico-diagnostico in vitro / Dispositif médical de diagnostic in vitro / Dispositivo médico para diagnóstico in vitro / In-Vitro-Diagnostikum	EC REP	Authorized Representative in the European Community / Rappresentante Autorizzato nella Comunità Europea / Mandataire dans la Communauté européenne / Representante autorizado en la Comunidad Europea / Bevollmächtigter in der Europäischen Gemeinschaft
CE	This product fulfills the requirements of Directive 98/79/EC on in vitro diagnostic medical devices / Questo prodotto soddisfa i requisiti della Direttiva 98/79/CE sui dispositivi medico-diagnostici in vitro / Ce produit répond aux exigences de la Directive 98/79 CE relative aux dispositifs médicaux de diagnostic in vitro / Este producto cumple con las exigencias de la Directiva 98/79/CE sobre los productos sanitarios para diagnóstico in vitro / Dieses Produkt entspricht den Anforderungen der Richtlinie über In Vitro Diagnostica 98/79/EG.	SMP PREP DIL SPE	Sample Preparation Apparatus containing Sample Diluent / Dispositivo per la preparazione del campione contenente il diluente del campione / Système pour la préparation de l'échantillon, diluant inclus / Aparato para Preparación de Muestra con Diluyente de Muestra / System zur Probenbereitung, in dem sich Probenverdünnungspuffer befindet
REF	Catalogue number / Numero di catalogo / Référence du catalogue / Número de catálogo / Bestellnummer		Do not freeze / Non congelare / Ne pas congeler / No congelar / Nicht Einfrieren
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Gebrauchsanweisung beachten	RoHS	Restriction of Hazardous Substances / Restrizione all'uso di sostanze pericolose / Limitation of substances dangereuses / Restricción de Substancias Nocivas / Beschränkung der Verwendung bestimmter gefährlicher Stoffe
	Manufacturer / Fabbricante / Fabricant / Fabricante / Hersteller		Caution - consult accompanying documents / Attention, vedere le istruzioni per l'uso / Attention voir notice d'instructions / Atención, ver instrucciones de uso / Achtung, Begleitdokumente beachten
	Contains sufficient for <n> tests / Contiene sufficiente per "n" saggi / Contenu suffisant pour "n" tests / Contenido suficiente para <n> ensayos / Inhalt ausreichend für <n> Prüfungen	BUF RXN	Reaction Buffer / Tampone di reazione / Solution de réaction tamponnée / Tampón de Reacción / Reaktionspuffer
	Temperature limitation / Limiti di temperatura / Limites de température / Limite de temperatura / Temperaturbegrenzung		ETL Registered Mark Certified / Marchio di certificazione registrato a livello nazionale / Certifíe Conforme ETL / Marca de Certificación Registrada Nacional / ETL Konform beglaubigt
SN	Serial number / Numero di serie / Numéro de série / Número de serie / Seriennummer		Recycle - do not dispose of as general waste / Reciclare - non eliminare come rifiuto generico / Recycler - ne pas jeter dans une poubelle / Recycle - no desecho como basura general / Recycling- dieses Produkt nicht über den Hausmüll entsorgen
TEST	Test Device / Dispositivo test / Dispositif de test / Dispositivo de Prueba / Testgerät	HT TUBE	Heat Treatment Tube / Provetta per il trattamento termico / Tube pour le traitement thermique / Tubo de tratamiento de calor / Röhren zur Hitzebearbeitung
	Date of manufacture / Data di fabbricazione / Date de fabrication / Fecha de fabricación / Herstellungsdatum		For MD Performance Evaluation Only / Soltanto per valutazione delle prestazioni / Réactifs IVD réservés à l'évaluation des performances / Solo para evaluación del funcionamiento / Nur zur IVD Leistungsbeurteilung
	LASER RADIATION: Avoid Exposure to Beam / RADIAZIONE LASER: Evitare l'esposizione al raggio / RAYONNEMENT LASER: Éviter toute exposition au faisceau / Radiación Laser: Evite Exposición a los Rayos / LASERSTRAHLUNG: Direkten Kontakt mit dem Strahl vermeiden		HOT SURFACE: Keep hands Away from Hot Surfaces / Superficie calda: tenere le mani lontane dalle superfici calde / SURFACES CHAUDES: Ne pas toucher les surfaces chaudes / Superficie Caliente: Mantenga las manos alejadas de la superficie caliente / Heiße Oberfläche: Kontakt mit heißen Oberflächen vermeiden
	CAUTION: Laser Radiation / ATTENZIONE: Radiazione Laser / AVERTISSEMENT: Rayonnement Laser / Precaución: Radiación Laser / WARNUNG: Laserstrahlung	IPX-0	CAUTION: Protect from water / ATTENZIONE: Proteggere dall'acqua / AVERTISSEMENT: Protéger de l'humidité / Precaución: Proteja del agua / WARNUNG: vor Feuchtigkeit schützen
	CAUTION: Risk of Danger / ATTENZIONE: Pericolo / AVERTISSEMENT: Risque de danger / Precaución: Peligroso / WARNUNG: Riskogefahr	CONTROL	Assay Control / Controllo del test / Test de contrôle / Control de Ensayo / Kontrolltest
BUF	Buffer / Soluzione tampone / Solution tamponnée / Tampón / Puffer	MIN OIL	Mineral Oil / Olio Minerale / Huile Minerale / Aceite Mineral / Mineralöl
CONJ	Conjugate / Coniugato / Conjugué / Conjugado / Konjugat		Warning / Avvertenze / Mise En Garde / Advertencia / Warnhinweise
ST TUBE	Screw Top Tube / Provetta con tappo a vite / Tube à bouchon vissé / Tubo con tapa de rosca / Röhren mit Schnappverschluss	DET REAG	Detection Reagent / Reagente Diretto / Réactif de Détection / Reactivo de Detección / Nachweis Reagenz

For technical assistance, call Technical Support Services at 800-343-3858 between the hours of 8AM and 6PM, USA Eastern Standard Time. To place an order, call Customer Service Department at 800-543-1980.