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Nom. : 2000030 - 09/08

## STAINING OF PARASITIC ELEMENTS IN STOOLS

### INTENDED USE:

**IODESINE-COLOR FUMOUCZE®** is a method of fixation and coloration in tube that allows a good observation of nuclear structures (chromatin - caryosome), necessary to identify many protozoa.

**IODESINE Coloration technique replaces the M.I.F. Coloration one while using the same staining agents.** The solution contains neither formalin nor mercury derivative so as to improve the work station safety.

### PRINCIPLE:

During the parasitologic examination of stools, detection and identification of amoeba and / or other protozoa are facilitated by vegetative and cystic form coloration.

**IODESINE-COLOR FUMOUCZE®** allows the parasite coloration and makes a slight concentration of parasitic elements at the surface of the sediment. It also allows to delay the observation of the direct examination and to preserve the parasitic elements for a longer time.

### KIT CONTENT:

- Base for IODESINE
- Lugol
- Stirrers
- Package insert

### MATERIAL REQUIRED BUT NOT PROVIDED:

- Micropipette 10-50 µL, pipette delivering 3 mL
- Haemolysis tubes
- Pasteur pipettes
- Slides + coverglasses for microscopy
- Microscope
- Vortex
- Container for contaminated wastes

### STORAGE CONDITIONS:

Reagents are ready-to-use.

Store at +18°...+25°C, **sheltered from sunlight**, until the expiry date indicated on the box. Do not freeze.

### SAMPLES STORAGE:

Due to the fragility of some of parasitic stages as protozoa vegetative forms, it is recommended to treat stools as soon as possible after their collection.

### WARNINGS AND PRECAUTIONS:

- For in vitro diagnostic use.

- Only for professional use.
- Follow the instructions for use.
- In case of accidental spill of reagents, clean the surface with absorbent paper and water. In case of spill of sample, clean the surface with absorbent paper and bleach.
- Avoid contact of reagent with skin, eyes and mucous membranes. Do not ingest.
- The samples, reagents as well as contaminated materials and products must be eliminated in a container for contaminated wastes, according to the prevailing recommendations and regulations.
- **Base for IODESINE:** R 10: Flammable.

### TEST PROCEDURE:

#### a. IODESINE solution preparation :

In a haemolysis tube, place, using a micropipette, 25 µL of lugol and 3 mL of Base for IODESINE.

- Homogenize the stools.
- Take out a volume of stools equivalent to the size of a pea, and place it in the previously prepared IODESINE solution.
- Triturate the stools by using a stirrer and shake the mixture vigorously until a homogenous suspension is obtained (Vortex shaker).
- Leave to settle 24 hours at room temperature, sheltered from the light.
- Collect 2 samples: the first one at the surface of the sediment and the second one after having put this one in suspension.
- Place 2 samples on a slide, place a coverglass over and examine under a microscope. When parasites colouring is too important, it is recommended to increase the white light intensity for a good microscopic observation of these last ones.

### INTERPRETATION OF THE RESULTS:

The vegetative and cystic forms of protozoa appear colored in pink brown, more or less dark. The observation of helminth eggs and larva is not disturbed by the IODESINE coloration.

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