

# illumipro-10™

Automated isothermal amplification and detection system that provides fast, accurate results in a flexible package that fits into any lab.



- Optimal flexibility for stat or batch runs
- True walk away capability for optimized workflow
- Very low maintenance requirements enhance efficiency
- Instrument verifications can be performed with in-house personnel

# illumipro-10™

**Digital display** guides users through process

**Automatic system self-check**

**Barcode reader** for sample ID entry

**Reaction blocks** operate independently


**Small footprint** of instrument is only 8.5" x 11"



## Procedure



1. Load samples



2. Press "RUN"



3. Enter patient ID

## Ordering information

<i>illumipro-10</i> ™ Incubator Reader.....	610172
<i>illumipro-10</i> ™ Printer.....	610173
<i>illumipro-10</i> ™ Keyboard.....	610174
<i>illumipro-10</i> ™ Operations Manual.....	11007

**i** For more information, contact an *illumigene*® specialist at 1-888-763-6769 or visit us on the web at [www.meridianbioscience.com](http://www.meridianbioscience.com).